

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
26 August 2004 (26.08.2004)

PCT

(10) International Publication Number
WO 2004/072230 A2

- (51) International Patent Classification⁷: **C12N**
- (21) International Application Number:
PCT/US2004/002012
- (22) International Filing Date: 10 February 2004 (10.02.2004)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data:
10/361,004 10 February 2003 (10.02.2003) US
- (71) Applicant (for all designated States except US): **CLEAR-ANT, INC.** [US/US]; 11111 Santa Monica Boulevard, Suite 650, Los Angeles, CA 90025 (US).
- (72) Inventors; and
- (75) Inventors/Applicants (for US only): **McKENNEY, Keith** [US/US]; 11918 Glen Mill Road, Potomac, MD 20854 (US). **GILLMEISTER, Lidja** [US/US]; 9419 Lee Highway, Fairfax, VA 22031 (US). **MARLOWE, Kristina** [US/US]; 9419 Lee Highway, Fairfax, VA 22031 (US). **ARMISTEAD, David** [US/US]; 1810 North Wayne Street, Arlington, VA 22201 (US).
- (74) Agents: **McPHAIL, Donald, R.** et al.; Fleshner & Kim, LLP, P.O. Box 221200, Chantilly, VA 20153-1200 (US).
- (81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NA, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.
- (84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).
- Published:
— without international search report and to be republished upon receipt of that report
- For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.



WO 2004/072230 A2

(54) Title: REAL-TIME POLYMERASE CHAIN REACTION USING LARGE TARGET AMPLICONS

(57) Abstract: The present invention relates to methods for analyzing a target nucleic acid sequence in a biological material. More particularly, the present invention relates to methods for analyzing a target nucleic acid sequence by real time polymerase chain reaction using nucleic acid primers that are separated by at least about 750 nucleic acid residues in the target sequence.

REAL-TIME POLYMERASE CHAIN REACTION USING LARGE TARGET AMPLICONS

BACKGROUND OF THE INVENTION

5 1. Field of the Invention

The present invention relates to methods for analyzing a target nucleic acid sequence in a biological material. More particularly, the present invention relates to methods for analyzing a target nucleic acid sequence by real time polymerase chain reaction using nucleic acid primers that are separated by at least about 750 nucleic acid residues in the target
10 sequence.

2. Background of the Related Art

PCR (polymerase chain reaction) is a method for increasing the concentration of a segment of a target sequence in a mixture of nucleic acid sequences without cloning or
15 purification. (*See* K. B. Mullis *et al.*, U.S. Pat. Nos. 4,683,195 and 4,683,202).

This process for amplifying the target sequence consists of introducing two oligonucleotide primers to the sample containing the desired target nucleic acid sequence, followed by thermal cycling in the presence of a DNA polymerase. The two primers are complementary to their respective strands of the target sequence. To effect amplification,
20 the genetic material within the sample is first denatured and then the primers are annealed to their complementary sequences within the target molecule. Following annealing, the primers are extended with a polymerase so as to form a new pair of complementary strands.

The steps of denaturation, annealing and extension can be repeated many times (*i.e.*, denaturation, annealing and extension constitute one "cycle"; there can be numerous

"cycles") to obtain a high concentration of an amplified segment of the desired target sequence. The length of the amplified segment of the desired target sequence is determined by the relative positions of the primers with respect to each other, and therefore, this length is a controllable parameter. Because the desired amplified segments of the target sequence
5 become the predominant sequences (in terms of concentration) in the mixture, they are said to be "PCR amplified".

With PCR, it is possible to amplify a single copy of a specific target sequence in genomic DNA to a level detectable by several different methodologies (*e.g.*, hybridization with a labelled probe; incorporation of biotinylated primers followed by avidin-enzyme
10 conjugate detection; incorporation of ^{32}P -labelled deoxynucleotide triphosphates, *e.g.*, dCTP or dATP, into the amplified segment). In addition to genomic DNA, any oligonucleotide sequence can be amplified with the appropriate set of primer molecules.

End-point PCR is a polynucleotide amplification protocol. The amplification factor that is observed is related to the number (n) of cycles that have occurred and the efficiency
15 of replication at each cycle (E), which, in turn, is a function of the priming and extension efficiencies during each cycle. Amplification has been observed to follow the form E^n , until high concentrations of the PCR product have been made.

At these high product concentrations, the efficiency of replication tends to drop significantly. It has been suggested that this is probably due to the displacement of the
20 primers by the longer complementary strands of the PCR product. At concentrations in excess of 10^{-8} M, the rate of the two complementary PCR amplified product strands finding each other during the priming reactions becomes sufficiently fast that it may occur before or

concomitantly with the extension step of the PCR process. This ultimately leads to a reduced priming efficiency, and, consequently, a reduced cycle efficiency. Continued cycles of PCR lead to declining increases of PCR product molecules, until the PCR product eventually reaches a plateau concentration (the "end-point"), usually a concentration of approximately 10^{-8} M. As a typical reaction volume is about 100 microliters, this corresponds to a yield of about 6×10^{11} double stranded product molecules.

Real-time PCR is also a polynucleotide amplification protocol, but PCR product analysis occurs simultaneously with amplification of the target sequence. Detecting agents, such as DNA dyes or fluorescent probes, can be added to the PCR mixture before amplification and used to analyze PCR products during amplification. Sample analysis occurs concurrently with amplification in the same tube within the same instrument. This combined approach decreases sample handling, saves time, and greatly reduces the risk of product contamination, as there is no need to remove the samples from their closed containers for further analysis. The concept of combining amplification with product analysis has become known as "real-time" or "quantitative" PCR. (*See, e.g.*, WO/9746707A2, WO/9746712A2 and WO/9746714A1).

Originally, monitoring fluorescence each cycle of PCR involved the use of ethidium bromide. *See* Higuchi *et al.*, "Simultaneous amplification and detection of specific DNA sequences," *Bio/Technology* **10**:413-417 (1992); Higuchi *et al.*, "Kinetic PCR analysis: real time monitoring of DNA amplification reactions," *Bio/Technology* **11**:1026-1030 (1993). In that system, fluorescence was measured once per cycle as a relative measure of product concentration. Ethidium bromide detects double stranded DNA; thus, if the desired target nucleic acid sequence is present, fluorescence intensity increases with temperature cycling

(otherwise no fluorescence). Furthermore, the cycle number where an increase in fluorescence is first detected increases inversely proportionally to the log of the initial target sequence concentration. Other fluorescent systems have since been developed that are capable of providing additional data concerning the nucleic acid concentration.

5 A significant limitation in the use of real-time PCR is the length of the target nucleic acid sequence. That is, as the target amplicon length increase, the efficiency of real-time PCR decreases. Practical limits for target amplicon length in most commercially available PCR systems are generally less than 500 bp, usually in the range of 80-200 bp. Larger amplicons have been obtained by some, but to date there remains a need for routinely
10 amplifying large target sequences in real time PCR.

Each of the above references is incorporated by reference herein where appropriate for teachings of additional or alternative details, features and/or technical background.

SUMMARY OF THE INVENTION

15 An object of the invention is to solve at least the problems and/or disadvantages of the relevant art, and to provide at least the advantages described hereinafter.

Accordingly, it is an object of the present invention to provide methods for analyzing a target nucleic acid sequence by real time polymerase chain reaction using nucleic acid primers that are separated by at least about 750 nucleic acid residues in the target sequence.
20 Other objects, features and advantages of the present invention will be set forth in the detailed description of preferred embodiments that follows, and in part will be apparent from the description or may be learned by practice of the invention. These objects and

advantages of the invention will be realized and attained by the compositions and methods particularly pointed out in the written description and claims hereof.

In accordance with these and other objects, a first embodiment of the present invention is directed to a method for analyzing a target nucleic acid sequence, comprising: (i) adding to a biological material an effective amount of at least two nucleic acid primers that hybridize under stringent conditions to predetermined sequences of the target sequence and are separated by at least about 750 nucleic acid residues, (ii) amplifying the target nucleic acid sequence by a polymerase chain reaction which comprises adding a polymerase to the biological material and primers to form an amplification mixture and thermally cycling the amplification mixture between at least one denaturation temperature and at least one elongation temperature, and (iii) detecting and quantifying said target nucleic acid sequence. According to this embodiment of the present invention, during the thermal cycling, the elongation temperature is not more than about 70°C and the denaturation temperature is not more than about 95°C, and the amplification mixture is maintained at the denaturation temperature for a period of not more than about 30 seconds and at the elongation temperature for a period of not less than about 1 minute.

Additional advantages, objects, and features of the invention will be set forth in part in the description which follows and in part will become apparent to those having ordinary skill in the art upon examination of the following or may be learned from practice of the invention. The objects and advantages of the invention may be realized and attained as particularly pointed out in the appended claims.

BRIEF DESCRIPTION OF THE DRAWINGS

The invention will be described in detail with reference to the following drawings in which like reference numerals refer to like elements wherein:

Figure 1 shows forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of human Parvovirus B19 (SEQ ID NO.: 1).

Figure 2 shows forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of hepatitis B virus (SEQ ID NO.: 2).

Figure 3 shows forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of porcine parvovirus (SEQ ID NO.: 3).

Figure 4 shows forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of Sindbis virus (SEQ ID NO.: 4).

Figure 5 shows forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of West Nile virus (SEQ ID NO.: 5).

Figures 6A and 6B show forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of the 16S ribosomal RNA gene (SEQ ID NO.: 6) and the 23S ribosomal RNA gene of *Escherichia coli* (SEQ ID NO.: 7).

Figures 7A and 7B show forward and reverse primers useful in preparing large target amplicons based on the genomic nucleic acid sequence of the 18S ribosomal RNA gene (SEQ ID NO.: 8) and the 25S ribosomal RNA gene of yeast (*S. cerevisiae*) (SEQ ID NO.: 9).

Figure 8 shows forward and reverse primers useful in preparing large target amplicons based on the nucleic acid sequence of human mitochondrial DNA (SEQ ID NO.: 10).

5 **DETAILED DESCRIPTION OF PREFERRED EMBODIMENTS**

A. *Definitions*

Unless defined otherwise, all technical and scientific terms used herein are intended to have the same meaning as is commonly understood by one of ordinary skill in the relevant art.

10 As used herein, the singular forms "a," "an," and "the" include the plural reference unless the context clearly dictates otherwise.

As used herein, the term "biological material" is intended to mean any substance derived or obtained from a living organism. Illustrative examples of biological materials include, but are not limited to, the following: cells; tissues; blood or blood components; 15 proteins, including recombinant and transgenic proteins, and proteinaceous materials; enzymes, including digestive enzymes, such as trypsin, chymotrypsin, alpha-galactosidase and iduronate-2-sulfatase; immunoglobulins, including mono and polyimmunoglobulins; botanicals; food and the like. Preferred examples of biological materials include, but are not limited to, the following: ligaments; tendons; nerves; bone, including demineralized bone 20 matrix, grafts, joints, femurs, femoral heads, etc.; teeth; skin grafts; bone marrow, including bone marrow cell suspensions, whole or processed; heart valves; cartilage; corneas; arteries and veins; organs, including organs for transplantation, such as hearts, livers, lungs, kidneys, intestines, pancreas, limbs and digits; lipids; carbohydrates; collagen, including native,

afibrillar, atelomeric, soluble and insoluble, recombinant and transgenic, both native sequence and modified; chitin and its derivatives, including NO-carboxy chitosan (NOCC); stem cells, islet of Langerhans cells and other cells for transplantation, including genetically altered cells; red blood cells; white blood cells, including monocytes; and platelets.

- 5 Additional examples of biological materials include forensic samples, human or animal remains, stomach contents, mummified remains of a once-living organism, fossilized remains, a product of manufacture containing or previously in contact with a biological material, and fomites.

10 **B. *Particularly Preferred Embodiments***

A first particular preferred embodiment of the present invention is directed to a method for analyzing a target nucleic acid sequence in a biological material, comprising:

- (i) adding to a biological material an effective amount of at least two nucleic acid primers, wherein these nucleic acid primers hybridize under stringent conditions to two
15 predetermined nucleic acid sequences of the target nucleic acid sequence that are separated by at least about 750 nucleic acid residues,
- (ii) amplifying the target nucleic acid sequence by a polymerase chain reaction which comprises adding a polymerase to the biological material and primers to form an amplification mixture and then thermally cycling the amplification mixture between at least
20 one denaturation temperature and at least one elongation temperature; and
- (iii) detecting and quantifying said target nucleic acid sequence.

According to this preferred embodiment of the present invention, the elongation temperature is not more than about 70°C and the denaturation temperature is not more than about 95°C. Additionally, according to this preferred embodiment of the present invention, during each thermal cycle, the amplification mixture is maintained at the denaturation temperature for a period of not more than about 30 seconds and at the elongation temperature for a period of not less than about 1 minute.

According to preferred embodiments of the present invention, the target nucleic acid sequence preferably contains between about 500 and about 50,000 nucleic acid residues. More preferably, the target nucleic acid sequence contains between about 1000 and about 10,000 nucleic acid residues, even more preferably between about 2000 and about 5000 nucleic acid residues and most preferably between about 2500 and about 5000 nucleic acid residues.

The nucleic acid primers are each selected based on their ability to generate the desired target nucleic acid sequence under the appropriate PCR conditions. Accordingly, each primer must be specific for the desired target nucleic acid sequence. Similarly, each primer must be selected so that they are not self-complementary or complementary to another primer (or probe, if present).

According to preferred embodiments of the present invention, the sequences on the target sequence that correspond to the two primer sequences are separated by at least 750 nucleic acid residues. Preferably, the sequences which correspond to the primers are separated by at least about 1000 nucleic acid residues, more preferably at least about 2000 nucleic acid residues, even more preferably at least about 3000 nucleic acid residues, still

even more preferably at least about 4000 nucleic acid residues and most preferably at least about 5000 nucleic acid residues. According to an alternative embodiment of the present invention, the sequences on the target sequence that correspond to the two primer sequences are separated by only about 500 nucleic acid residues.

5 The polymerase chain reaction employed in the inventive methods is performed according to the methods and techniques known to those skilled in the art, *i.e.*, a nucleic acid primer pair is added to the biological material containing the sequence of interest to form an amplification mixture that is then thermally cycled for a sufficient period of time to amplify the desired sequence. The thermal cycling generally comprises cycling the amplification
10 mixture between at least one denaturation temperature and at least one elongation temperature. Preferably, the thermal cycling comprises cycling the amplification mixture between at least one denaturation temperature, at least one annealing temperature and at least one elongation temperature.

Specific temperatures for use in denaturation, elongation and/or annealing may be
15 determined empirically by one skilled in the art based, for example, on the specific target sequence being amplified and the particular probes employed. Likewise, the specific time(s) that the amplification mixture is maintained at the various denaturation, elongation and/or annealing temperature(s) may be determined empirically by one skilled in the art based on similar considerations.

20 According to particularly preferred embodiments of the present invention, the elongation temperature selected for use in the PCR of the inventive methods is not more than about 70°C. More preferably, the elongation temperature selected is between about

60°C and about 69°C, and even more preferably between about 65°C and about 69°C. Most preferably, the elongation temperature employed in the PCR of the inventive methods is about 68°C.

According to additional preferred embodiments of the present invention, the
5 denaturation temperature selected for use in the PCR of the inventive methods is not more than about 95°C. More preferably, the denaturation temperature selected is between about 90°C and about 95°C, and even more preferably between about 93°C and about 95°C. Most preferably, the denaturation temperature employed in the PCR of the inventive methods is about 94°C.

10 According to other preferred embodiments of the present invention, when the thermal cycling includes an annealing temperature, the annealing temperature selected is about 5-10°C below the melting temperature of the primers being employed. Preferably, the annealing temperature selected is not more than about 65°C. More preferably, the annealing temperature selected is between about 57°C and about 63°C, and even more preferably
15 between about 58°C and about 62°C. Most preferably, the annealing temperature employed in the PCR of the inventive methods is about 60°C.

According to additional preferred embodiments of the present invention, during each thermal cycle, the amplification mixture is maintained at the elongation temperature for a period of not less than about 1 minute. More preferably, during each thermal cycle, the
20 amplification mixture is maintained at the elongation temperature for a period of not less

than about 2 minutes, and even more preferably for a period of not less than about 3 minutes.

According to particularly preferred embodiments of the present invention, the amplification mixture is maintained at the elongation temperature for a period of not less than about 2 minutes during the first cycle of the thermal cycling, and then the period during which said amplification mixture is maintained at the elongation temperature is increased by a period of about 5 seconds for each successive thermal cycle. Thus, for example, according to such embodiments of the present invention, if the amplification mixture was maintained at the elongation temperature for a period of 2 minutes during the first cycle of the thermal cycling, it would be maintained at the elongation temperature for a period of 2 minutes, 5 seconds for the second cycle, 2 minutes, 10 seconds for the third cycle, 2 minutes, 15 seconds for the fourth cycle, and so on until the thermal cycling is completed.

According to additional preferred embodiments of the present invention, during each thermal cycle, the amplification mixture is maintained at the denaturation temperature for a period of not more than about 1 minute. More preferably, during each thermal cycle, the amplification mixture is maintained at the denaturation temperature for a period of not more than about 45 seconds, and even more preferably for a period of not more than about 30 seconds, and still even more preferably for a period of not more than about 20 seconds. Most preferably, during each thermal cycle, the amplification mixture is maintained at the denaturation temperature for a period of not more than about 15 seconds, such as a period of about 10 seconds.

According to still other preferred embodiments of the present invention, when the thermal cycling includes an annealing temperature, the amplification mixture is maintained at the annealing temperature for a period of not less than about 30 seconds. More preferably, according to such embodiments, during each thermal cycle, the amplification mixture is
5 maintained at the annealing temperature for a period between 30 seconds and 2 minutes, and even more preferably for a period of not less than about 45 seconds. Most preferably, during each thermal cycle, the amplification mixture is maintained at the annealing temperature for a period of about 1 minute.

The number of thermal cycles employed in the PCR of the inventive methods may be
10 determined empirically by one skilled in the art depending, for example, on the suspected concentration of the target sequence of interest in the biological material being tested. According to preferred embodiments of the present invention, the amplification mixture is subjected to at least about 30 cycles of thermal cycling, and even more preferably at least about 40 cycles. Most preferably, the amplification mixture is subjected to at least about 50
15 cycles of thermal cycling.

The polymerase employed in the PCR of the inventive methods may be any of the suitable polymerases known to those skilled in the art. Preferably, the polymerase employed is a thermostable polymerase, *i.e.* a polymerase that is not adversely affected by the higher temperatures involved in thermal cycling. More preferably, the polymerase may be a *Taq*
20 polymerase, or a suitable derivative thereof and/or a proof-reading polymerase.

According to particularly preferred embodiments of the present invention, at least two polymerases are employed in the PCR of the inventive methods. Preferably, at least one

of the polymerases is a *Taq* polymerase or a suitable derivative thereof, such as TaqMan DNA polymerase (available from Applied BioSystems), and the other polymerase is a proof-reading polymerase, such as ProofStart DNA polymerase (available from Qiagen).

According to certain preferred embodiments of the present invention, the
5 amplification mixture further contains at least one thermostable inorganic pyrophosphatase. Suitable amounts of thermostable inorganic pyrophosphatase may be determined empirically by one skilled in art. Generally, when present, the ratio of thermostable inorganic pyrophosphatase to *Taq* polymerase is at least about 1:20, more preferably at least about 1:10 and even more preferably at least about 1:5.

10 The remaining parameters employed in the PCR of the inventive methods, such as the primer concentration (generally about 100-500 nM and preferably about 200 nM), magnesium concentration (generally 1.5-6 mM and preferably about 1.5 mM of magnesium sulfate and/or magnesium chloride), deoxyribonucleotide triphosphates (dNTP)
concentration (generally about 0.2-0.4 mM each and preferably about 0.2 mM each), probe
15 concentration (if present, generally about 50-800 nM, and preferably about 100 nM), may each be determined empirically by one skilled in the art using any of the known methods and techniques.

According to certain particularly preferred embodiments of the present invention, the deoxyribonucleotide triphosphates (dNTP) that are employed in the PCR of the inventive
20 methods are selected from the group consisting of C, T, G and A. Preferably, substantially no dUTP is present in the amplification mixture of the inventive methods. According to still

further preferred embodiments, substantially no uracil N-glycosylase is present in the amplification mixture of the inventive methods.

According to certain particularly preferred embodiments of the present invention, the amplification mixture further comprises at least one buffer solution. Suitable buffer solutions
5 are known and available to those skilled in the art. Particularly preferred buffer solutions include pH modifying buffers, such as buffers containing Tris-HCl, and buffers which maintain salt concentration, particular magnesium concentration, such as buffers containing KCl and/or $(\text{NH}_4)_2\text{SO}_4$.

After amplification using PCR, the first and second target nucleic acid sequences are
10 detected and quantified. This detecting and quantifying may be conducted using any of the methods and techniques known to those skilled in the art. For example, detecting and quantifying of the first and second nucleic acid sequences may be conducted by adding a suitable detecting agent, such as an intercalating dye, directly to the amplification mixture or by adding a suitable nucleic acid probe to the mixture, preferably either a suitable nucleic
15 acid probe in combination with a detecting agent or a suitable nucleic acid probe having a detectable label covalently or ionically attached thereto or complexed therewith.

Preferably, the target nucleic acid sequence is detected by adding at least one nucleic acid probe to the biological material being tested. Any nucleic acid probe employed in the inventive methods should contain sufficient nucleic acid residues to hybridizes selectively
20 under stringent conditions to a specific desired nucleic acid sequence, *i.e.* suitable probes will generally contain at least 16 nucleic acid residues, and preferably hybridizes selectively under stringent conditions to a specific nucleic acid sequence of the target nucleic acid sequence

that is not the same as the nucleic acid sequence of any of the primers. Suitable nucleic acid probes include, but are not limited to, 5' nuclease probes, hairpin probes, adjacent probes, sunrise probes and scorpion probes.

EXAMPLES

The following examples are illustrative, but not limiting, of the present invention. Other suitable modifications and adaptations are of the variety normally encountered by those skilled in the art and are fully within the spirit and scope of the present invention.

5 Example 1

Purpose: To demonstrate linear amplification of B19 DNA.

Materials: 1. B19 virus, titer 7.6×10^{11} iu/ml from Bayer;

2. SNAP whole blood DNA isolation kit;

3. Forward Primer: Prism 5 (Figure 1) (SEQ ID NO.: 18);

10 4. Reverse Primer: Prism 6 (Figure 1) (SEQ ID NO.: 20);

5. Probe 3 (Figure 1) (SEQ ID NO.: 19) labeled with FAM at 5'

end and TAMRA at 3' end;

6. TaqMan Universal Master Mix, (ABI; cat. no. 4304437);

7. DNASE, RNASE free water;

15 8. ABI 96 well plate and adhesive cores;

9. ANI 7000.

Procedure: 1. Followed SNAP protocol for extraction of 100 μ l B19

sample, eluted in 100 μ l TE;

2. Diluted primers to 18 μ M with TE;

20 3. Diluted probe to 5 μ M with TE;

4. Prepared the following master mix:

TaqMan Master Mix: 25 µl;
Prism 5 (SEQ ID NO.: 18) 2.5 µl;
Prism 6 (SEQ ID NO.: 20) 2.5 µl;
Taqman Probe 2.5 µl;
5 Water: 12.54 µl;

5. Added 45 µl of master mix per well;
6. Serially diluted B19 DNA, adding water to the NTC well;
7. Sealed and centrifuged the plate at 2300 rpm for about 30

seconds;

- 10 8. Ran PCR program for 50 cycles.

Results: A standard dilution curve was observed for B19 infected plasma, validating primer pair Prism 5 and Prism 6 (SEQ ID NOS.: 18 and 20) with Probe 3 (SEQ ID NO.: 19).

15 Example 2

Purpose: To examine irradiated and unirradiated samples containing PPV using a 549 bp amplicon.

Materials: 1. PPV (irradiated at 0 kGy, 50 kGy, 65 kGy, 75 kGy or 85 kGy);

- 20 2. SNAP Protein Degradar;
3. Cell Lysis Buffer;

4. Tris-HCl;

5. Primers: Prism 11 and Prism 12 (Figure 3) (SEQ ID NOS.: 40 and 42, respectively); and

6. Probe 6 (Figure 3) (SEQ ID NO.: 41).

- 5 Procedure:
1. To 100 µl viral sample, added 50 µl tris-HCl buffer, 60 µl protein degrader, and 200 µl cell lysis buffer;
 2. Mixed and incubated for 25 minutes (5 minutes at 70°C);
 3. Diluted samples to 1/50, 1/500, 1/5000, 1/25000, 1/50000, 1/250000 and 1/500000;
 - 10 4. Ran PCR for 55 cycles.

Results: Results showed that unirradiated material had regular dilution series curves, irradiated material (50 kGy) behaved differently, dilute material did not amplify showing a reduction in the number of copies of the target sequence.

15 Example 3

Purpose: To determine effects of gamma irradiation (0 kGy sample, 50 kGy sample, mixture of 0+50kGy sample and 75 kGy sample) on samples containing PPV analyzed by PCR.

- Materials:
1. PPV (irradiated at 0 kGy, 50 kGy or 75 kGy);
 - 20 2. Primers: Prism 11 & Prism 12, Probe 6 (Figure 3) (SEQ ID NOS.: 40, 42 and 41, respectively);

3. Primers: Prism 1 & Prism 2, Probe 1 (Figure 3) (SEQ ID NOS.: 43, 45, and 44, respectively).

Procedure: 1. Diluted samples containing PPV to 1/100, 1/1000, 1-2000, 1/10000, 1/20000, 1/40000 and 1/400000 (0 kGy, 50 kGy, 0+50 kGy and 75 kGy);

5 2. Ran PCR program for 55 cycles.

Results: Irradiation to 50 kGy of PPV material reduced amplification of 549 bp amplicon.

Example 4

10 Purpose: To examine the relative effectiveness of Qiagen and Taqman reagents on samples containing PPV.

Materials: 1. PPV DNA (phenol extracted);
2. Taq PCR Core Kit;
3. ProofStart DNA polymerase;
15 4. Taqman Universal PCT Master Mix;
5. Prism 1, 2, 11 and 17 (Figure 3) (SEQ ID NOS.: 43, 45, 40, and 47 respectively);
6. Probes 1 and 6 (Figure 3) (SEQ ID NOS.: 44 and 41, respectively);

20 7. Agarose;
8. TAE;
9. EtBr.

Procedure:

1. Prepared the following four master mixes:

| | | |
|-------------------|----------------|---------------|
| a. Qiagen: | 1 | 2 |
| 10x buffer: | 30 μ l | 25 μ l |
| dNTP's: | 9 μ l | 7.5 μ l |
| pA: | 8.34 μ l | 6.95 μ l |
| pB: | 8.34 μ l | 6.95 μ l |
| taq: | 6 μ l | 5 μ l |
| H ₂ O: | 187.32 μ l | 156.1 μ l |
| probe: | 15 μ l | 12.5 μ l |
| b. Taqman: | 3 | 4 |
| Master Mix: | 150 μ l | 125 μ l |
| pA: | 15 μ l | 12.5 μ l |
| pB: | 15 μ l | 12.5 μ l |
| probe: | 15 μ l | 12.5 μ l |
| H ₂ O: | 69 μ l | 57.5 μ l |

2. Pipetted 44 μ l of master mix 1 into row D, wells 1 and 2;
row E, wells 1 and 2; and row H, well 1, of a well plate;

3. Pipetted 44 μ l of master mix 2 into row D, wells 3 and 4;
and row E, wells 3 and 4, of a well plate;

4. Pipetted 44 μ l of master mix 3, into row F, wells 1 and 2;
row G, wells 1 and 2; and row H, well 3, of a well plate;

5. Pipetted 44 µl of master mix 4 into row F, wells 3 and 4; and row G, wells 3 and 4, of a well plate;

6. Added 1 µl of ProofStart taq to row D, wells 1-4 and row F, wells 1-4 and added 1 µl water to remaining wells;

5 7. Added 5 µl water to row H, wells 1 and 3 and added 5 µl PPV DNA to remaining wells;

8. Ran PCR for 40 cycles.

Results: Qiaqen Master with ProofStart taq produced functional large amplicons in realtime PCR with PPV DNA more efficiently than the TaqMan master mix.

10

Example 5

Purpose: To examine the effects of proofstart in amplifying large amplicons and to examine the effects of 50 kGy irradiation on PPV.

Materials:

1. PPV DNA (irradiated to 0 kGy and 50 kGy);
- 15 2. Taq PCR Core Kit;
3. Proofstart DNA polymerase;
4. Prism 11, 16 and 17 (Figure 3) (SEQ ID NOS.: 40, 46 and 47, respectively);
5. Agarose;
- 20 6. Ethidium Bromide;
7. TAE buffer.

Procedure: 1. Set up PCR master mix as follows:

10x buffer: 50 µl
dNTP's: 15 µl
pA: 13.9 µl (primer 11)
taq: 10 µl
5 water: 347.2 µl

2. Placed aliquots into PCR tubes;
3. Added either primer 16 or 17 to PCR tubes;
4. Added PPV DNA (diluted to 1:100) to each PCR tube:
5. Added 10 µl proofstart to half of the samples (2 at 0 kGy
10 and 2 at 50 kGy);
6. Performed PCR (about 55 cycles)
7. Poured a 1% gel and ran at 100 V for 20 minutes.

Results: Addition of a proofreading polymerase resulted in improved
amplication of longer amplicons. Delay in amplification of target sequence in irradiated
15 samples is proportional to damage done to viral genetic material.

Example 6

Purpose: To examine the effect of TSP concentration on amplification of large
target amplicons in gamma irradiated and unirradiated PPV.

- 20 Materials:
1. TSP (cat. no. M02965);
 2. Qiagen Core kit;

3. ProofStart DNA polymerase;
4. PPV (irradiated to 0 kGy or 50 kGy).

Procedure:

1. Prepared a master mix (standard PCR set-up) for each (TSP Taq 1:20, 1:10, 1:5);
- 5 2. Added 43.61 µl of each master mix (TSP titration) to PCR tubes;
3. Added 1.39 µl of primers 16, 17 or 19 (Figure 3) (SEQ ID NOS.: 46, 47, and 49, respectively) to appropriate PCT tubes;
4. Added 5 µl water to the negative control, which
- 10 contained primer pair 11, 16 (Figure 3) (SEQ ID NOS.: 40 and 46, respectively);
5. Diluted PPV 1:100;
6. Added PPV to PCR tubes;
7. Performed PCR;
8. Poured a 1% gel and ran at 100 V for 20 minutes.

- 15 Results: Under these conditions, addition of TSP resulted in increased amplification of target amplicons in both irradiated and unirradiated samples, but irradiation of PPV resulted in decreased amplification of target amplicons.

Example 7

- 20 Purpose: To examine the effects of gamma irradiation on amplification of PPV target amplicons of various sizes.

Materials:

1. PPV DNA (irradiated to 0 kGy or 50 kGy);

2. Taq PCR Core Kit;
3. ProofStart DNA Polymerase;
4. Prism 11, 16, 17, 18 and 19 (Figure 3) (SEQ ID NOS.: 40, 46, 47, 48, and 49, respectively);

5. Agarose;
6. TAE;
7. Ethidium Bromide.

Procedure:

1. Prepared PCR Master Mix as follows:

| | | |
|----|------------|---------------|
| | 10x Buffer | 5 μ l |
| 10 | dNTPs | 1.5 μ l |
| | pA | 1.39 μ l |
| | pB | 1.39 μ l |
| | taq | 1 μ l |
| | water | 33.72 μ l |
| 15 | PPV | 5 μ l. |

2. Aliquoted samples into PCR tubes;
3. Ran PCR;
4. Poured a 1% agarose gel and ran at 120 V for about 1.5 hours.

20 Results: Irradiation to 50 kGy resulted in decreased amplification of larger target amplicons.

Example 8

Purpose: To examine PCR sensitivity and determine log reduction of PPV in samples irradiated to 50 kGy and having a starting concentration of 2.5×10^7 gEq.

Materials:

1. Standard PCR reagents (Qiacen Core Kit, TSP, Proofstart, etc.);
2. Primers 11 and 17 (Figure 3) (SEQ ID NOS.: 40 and 47, respectively);
3. PPV extract.

Procedure:

1. Prepared master mix with primers 11 and 17;
2. Performed a 10 fold dilution series from 10^7 to 10^0 of PPV extract;
3. Pipetted 45 μ l of master mix into PCR tubes;
4. Pipetted 5 μ l of each PPV dilution into appropriated PCR tubes;
5. Added 5 μ l water to control;
6. Ran PCR;
7. Ran samples in 1% agarose at 100V for about 47 minutes.

Results: Irradiation of sample to 50 kGy resulted in decreased amplification of target amplicon across all concentration ranges.

Example 9

Purpose: To examine PCR sensitivity and determine log reduction of PPV irradiated to 50 kGy and having a starting concentration of 2.5×10^7 gEq.

- Materials:
1. TSP;
 2. Standard PCR kit (Qiacen with ProofStart Polymerase);
 3. Primers 11 and 19 (Figure 3) (SEQ ID NOS.: 40 and 49, respectively);
 4. PPV Extract (Irradiated to 0 kGy and 50 kGy).

- Procedure:
1. Prepared master mix with primers 11 and 19 (SEQ ID NOS.: 40 and 49, respectively);
 2. Performed a 10 fold dilution series from 10^7 to 10^0 of PPV extract;
 3. Pipetted 45 μ l of master mix into PCR tubes;
 4. Pipetted 5 μ l of each PPV dilution into appropriate PCR tubes;
 5. Added 5 μ l water to control;
 6. Ran PCR as follows:
 - 95°C for 2 minutes (1 cycles)
 - 94°C for 10 seconds (40 cycles)
 - 60°C for 1 minute (40 cycles)
 - 68°C for 2 minutes (40 cycles);
 7. Cooled to 4°C;

8. Ran samples on 1% agarose gel in 1x TAE and 5 µl/100 ml ethidium bromide at 100 V for 52 minutes (5 µl on gel).

Results: Irradiation to 50 kGy resulted in decreased amplification of target amplicon across all concentration ranges. For unirradiated samples, relative band strength of
5 observed target amplicon decreased with decreasing concentration.

Example 10

Purpose: Primer validation for B19 using probe 7 (SEQ ID NO.: 12) and various primers.

10 Materials:

1. B19 IGIV Paste (irradiated to 0 kGY or 50 kGy);
2. EXB;
3. Proteinase;
4. yeast tRNA
5. phenol chloroform isoamyl alcohol;
- 15 6. 3M NaAc;
7. isopropanol;
8. 70% EtOH;
9. TE buffer;
10. Prisms 5, 6, 20, 21, 22, 23, 24, 25, 26 (Figure 1) (SEQ ID
- 20 NOS.: 18, 20, 11, 13, 14, 15, 16, 17, and 21, respectively);
11. Qiagen reagents;
12. Ampligold Taq;

13. ProofStart Polymerase;

14. Agarose;

15. TAE;

16. Ethidium Bromide.

5 Procedure:

1. Prepared a Master Mix as follows:

Buffer 5 µl

DNTP 1.5 µl

Taq 1 µl

DNA 5 µl

10 water 34.72 µl

2. Pipetted Master Mix into PCR tubes;

3. Added the following primer pairs to appropriate PCR tubes:

20&21 (SEQ ID NOS.: 11 and 13, respectively); 20&22 (SEQ ID NOS.: 11 and 14,
respectively); 20&23 (SEQ ID NOS.: 11 and 15, respectively); 20&24 (SEQ ID NOS.: 11
15 and 16, respectively); 20&25 (SEQ ID NOS.: 11 and 17, respectively); 20&6 (SEQ ID NOS.:
11 and 20, respectively); 20&26 (SEQ ID NOS.: 11 and 21, respectively); 5&6 (SEQ ID
NOS.: 18 and 20, respectively);

4. Ran PCR;

5. ran 1% gel for about 1 hour.

20 Results: All tested primers yielded desired target amplicons.

Example 11

Purpose: Use of PCR multiplexing with target amplicons of about 112 bp and about 2.4 kbp for B19 virus in samples irradiated to 0 kGy or 50 kGy.

- Materials:
1. TSP thermostable inorganic pyrophosphatase
 - 5 2. Standard PCR reagents;
 3. B19 viral extract (irradiated to 0 kGy and 50 kGy);
 4. Primers 5, 6, 20 and 25 (Figure 1) (SEQ ID NOS.: 18, 20, 11 and 17, respectively);
 5. Taq;
 - 10 6. ProofStart Polymerase.

- Procedure:
1. Prepared standard PCR set-up with 3x master mixes, for each primer set (primer sets: 5&6 (SEQ ID NOS.: 18 and 20, respectively); 20&25 (SEQ ID NOS.: 11 and 17, respectively); 5&6 (SEQ ID NOS.: 18 and 20, respectively); and 20&25 (SEQ ID NOS.: 11 and 17, respectively));
 - 15 2. Prepared appropriate PCR tubes containing the following primer pairs: (5, 6) 0 kGy; (5, 6) 50 kGy; (20, 25) 0 kGy; (20, 25) 50 kGy; (5, 6) & (20, 25), 0 kGy; and (5, 6) and (20, 25), 50 kGy;
 3. Added 5 µl B19 to PCR tubes containing 45 µl of appropriate master mix;
 - 20 4. Added 5 µl water to control;
 5. Ran PCR.

6. Ran samples on 1% agarose gel at 100 V for about 17 minutes.

Results: PCR multiplexing is effective for mixtures containing large target amplicons and small target amplicons. Irradiation to 50 kGy resulted in decreased
5 amplification of the large target amplicon relative to the small target amplicon.

Example 12

Purpose: Irradiated and unirradiated samples containing B19 viral material were examined using real time PCR.

10 Materials:

1. B19 viral material (irradiated to 0 kGy and 50 kGy);
2. Prism pairs (20, 21) (SEQ ID NOS.: 11 and 13, respectively)
and (20, 26) (SEQ ID NOS.: 11 and 21) (Figure 1);
3. Qiagen PCR reagents;
4. Qiagen ProofStart;
- 15 5. Agarose;
6. TAE (1x);
7. sample loading buffer (SLB).

Procedure:

1. Prepared standard samples containing primer pairs with 10^{11}
to 10^1 dilution series;
- 20 2. Ran PCR (40 cycles);
3. Ran gel on 1% agarose (8 μ l PCR product, 1 μ l SLB) at 100
V for about 20 minutes.

Results: Irradiation to 50 kGy resulted in decreased amplification of large target amplicon. Unirradiated samples exhibited a regular dilution pattern.

Example 13

5 Purpose: To investigate the effect of gamma irradiation on samples containing HBV and irradiated to 50 kGy.

Materials:

1. HBV (irradiated to 0 kGy and 50 kGy);
2. Taq PCR Core Kit;
3. ProofStart DNA polymerase;
- 10 4. Prisms 34, 9, 10, 15, 29, 30, 31, 36 and 37 (Figure 2) SEQ ID NOS.: 31, 22, 24, 25, 27, 32, 34, 28, and 29, respectively);
5. Agarose;
6. TAE Buffer;
7. ethidium bromide.

15 Procedure: 1. Prepared PCR master mix as follows:

| | |
|----------------------|---------------|
| 10x PCR buffer | 5 μ l |
| dNTPs | 1.39 μ l |
| primers | 1.39 μ l |
| taq | 1 μ l |
| 20 ProofStart | 1 μ l |
| water | 33.22 μ l |

TSP 0.5 µl

2. Aliquoted 43.61 µl of master mix into PCR tubes.

Appropriate tubes contained the following primer pairs: (3, 4); (9, 10); (9, 15); (9, 29); (9, 30); (9, 31); (36, 37); and (9, 31), for both irradiated and unirradiated samples;

- 5 3. Added 5 µl HBV per tube (irradiated or unirradiated);

4. Ran PCR as follows:

50°C for 2 minutes (one cycle)

95°C for 2 minutes (one cycle)

94°C for 10 seconds (40 cycles)

10 60°C for 1 minute (40 cycles)

68°C for 2 minutes, five seconds (40 cycles);

5. Ran 1% agarose gel (9 µl sample + 1 µl sample buffer) at

100v for about 20 minutes.

Results: Irradiated samples showed no amplification of large target amplicons,
15 indicating degradation of HBV genetic material by irradiation to 50 kGy.

Example 14

Purpose: To investigate the effect of gamma irradiation on samples containing
HBV DNA and irradiated to 50 kGy.

20 Materials:

1. HBV DNA material (irradiated to 0 kGy and 50 kGy);
2. Taq PCR Core Kit (Qiagen, cat. no. 201223);
3. ProofStart Taq Polymerase (Qiagen, cat. no. 20);

4. Prisms 10, 13, 30, 36 and 37 (Figure 2) (SEQ ID NOS.: 24, 26, 32, 28, and 29, respectively);

5. Agarose;

6. TAE Buffer;

5

7. Ethidium Bromide.

Procedure:

1. Prepared the following master mix:

| | | |
|----|----------------------------|-----------------|
| | 10x buffer | 60 μ l |
| | dNTP | 18 μ l |
| | primer 36 (SEQ ID NO.: 28) | 16.68 μ l |
| 10 | Taq | 12 μ l |
| | ProofStart | 12 μ l |
| | water | 440.64 μ l; |

2. Pipetted 46.61 μ l of master mix into PCR tubes;

3. Added 1.39 μ l of reverse primer (10, 13, 30 or 37) (SEQ ID NOS.: 24, 26, 32, and 29, respectively) and 2 μ l HBV DNA (0 kGy and 50 kGy) to appropriate tubes;

4. Ran PCR for 50 cycles;

5. Poured a 1% agarose gel (8 μ l PCR product + 1 μ l sample buffer) at 100 V for about 20 minutes.

20 Results: Irradiated samples showed no amplification of large target amplicons, indicating degradation of HBV DNA by irradiation to 50 kGy.

Example 15

Purpose: HBV amplification of nested primer set (about 80 bp, 400 bp and 697 bp) in samples containing ascorbate, including digestion of 0 kGy and 50 kGy samples with
5 exonuclease I prior to PCR amplification.

Materials:

1. HBV DNA (irradiated to 0 kGy and 50 kGy, with and without ascorbate);
2. Primer sets: (9, 10) (SEQ ID NOS.: 22 and 24, respectively); (9, 15) (SEQ ID NOS.: 22 and 25, respectively); and (9, 13) (SEQ ID NOS.: 22 and 26,
0 respectively) (Figure 2);
3. Exonuclease I;
4. Standard PCR reagents.

Procedure:

1. Diluted HBV samples to 1/500, 1/2000 and 1/10000;
2. Digested 1 µl raw HBV extract in 0.25 µl Exonuclease I, 10
5 µl 10x Exonuclease I buffer and 88.75 µl water at 37°C for 30 minutes, inactivated at 80°C
for 20 minutes;
3. Dilutes digested HBV to 1/2000 and 1/10000;
4. Ran 55 cycles PCR.

Results: Irradiated showed no amplification of large target amplicon (697 bp),
10 indicating degradation of HBV DNA by irradiation to 50 kGy.

Example 16

Purpose: To investigate the amount of bacterial and fungal DNA present in pulverized tendon samples.

- Materials:
1. E. Coli samples (tendon) – 0 or 50 kGy + stabilizer
5 (6.65x10¹⁰ CFU/μl);
 2. C. Albicans samples (tendon) – 0 or 50 kGy + stabilizer
(3.55x10⁹ CFU/μl);
 3. Staph. Aureus samples;
 4. Control tendon;
 5. Dneasey tissue kit (Qiagen, cat. no. 69504);
 6. Taq PCR Core Kit (Qiagen, cat. no. 201223);
 7. ProofStart Taq Polymerase (Qiagen, cat. no. 202205);
 8. Primers: Ribo 7 and 8, and Ribo 10, 11, 12, 13, 14 (Figures
6A and 6B) SEQ D NOS.: 69, 70, 71, 72, and 73, respectively); and Fungi 1, 2, 3, 4, 5, 6, 7, 8
5 (Figures 7A and 7B) (SEQ ID NOS.: 75, 77, 78, 79, 80, 81, 82, and 83 respectively);
 9. Probes: FAM-RIBO
Fungi Probe (Figure 7A) (SEQ ID NO.: 76) labeled with
FAM at 5' end and TAMRA at 3' end;
 10. Microcon YM Centrifugal Filter Unit;
- Procedure:
1. Using 0.05 tendon samples for E. coli and C. albicans,
followed the Qiagen extraction profile;
 2. Prepared the following master mixes:

| | | Mix 1 | Mix 2 |
|---|-------------------------|---------------|----------------|
| | 10x buffer | 150 μ l | 85 μ l |
| | dNTPs | 45 μ l | 25.5 μ l |
| | Ribo 7 | 41.7 μ l | --- |
| 5 | Fungi 1 (SEQ ID NO. 75) | --- | 23.65 μ l |
| | Taq | 30 μ l | 17 μ l |
| | ProofStart | 30 μ l | 17 μ l |
| | Water | 936.6 μ l | 530.74 μ l |
| 0 | FAM-RIBO | 75 μ l | --- |
| | Fungi Probe | --- | 42.5 μ l |

3. Filtered master mixes using Microcon filter units for 30

minutes at 100x g;

4. Pipetted 43.6 μ l of Mix 1 into: rows A-D, columns 1-6; rows
 5 A-C, column 9; and row E, column 12;

5. Pipetted 43.6 μ l of Mix 2 into: rows E-F, columns 1-7 and
 row H, column 12;

6. Pipetted 1.39 μ l of reverse primer into appropriate well;

7. Pipetted 5 μ l DNA into appropriate wells;

20 8. Ran PCR.

Results: Irradiation with 50 kGy resulted in decreased amplification of large target amplicons, indicating degradation of the pathogen genetic material caused by irradiation.

5 Example 17

Purpose: To show functionality of E. coli primers for RT-PCR using large target amplicons.

- Materials:
1. E. coli prepared from overnight culture;
 2. Dneasy Tissue Kit (Qiagen, cat. no. 96504);
 - 0 3. Taq PCR Core Kit (Qiagen, cat. no. 201223)
 4. ProofStart DNA polymerase (Qiagen, cat. no. 202205);
 5. Microcon YM-100 Centrifugal Filter Unit (cat. no. 42413);
 6. Primers: Ribo 1-6 (SEQ ID NOS.: 62, 64, 65, 66, 67, and
68, respectively) and Ribo 7-9;
 - 5 7. Agarose;
 8. TAE Buffer;
 9. Ethidium Bromide.

- Procedure:
1. Pipetted 1 ml of E. coli culture into each of 10 1.5 tubes;
 2. Centrifuged all 10 tubes for 5 minutes at maximum speed;
 - 0 3. Discarded supernatant;
 4. Placed 8 tubes in -80°C and used 2 tubes for extraction
following the Qiagen protocol;

5. Prepared Master Mix as follows:

| | |
|------------|--|
| 10x Buffer | 5 μ l |
| dNTPs | 1.5 μ l |
| pA | 1.39 μ l (Ribo 1 (SEQ ID NO.: 62)) |

5 or (Ribo 7)

pB 1.39 μ l (Ribo 2, 3, 4, 5, or 6) (SEQ ID

NOS.: 64, 65, 66, 67, or 68, respectively) or (Ribo 8 or 9)

Taq 1 μ l

ProofStart 1 μ l

0 Water 33.22 μ l

TSP 0.5 μ l

6. Mixed Master Mix by inversion;

7. Pipetted Master mix into a Microcon Centrifugal Filter Unit
and centrifuged for 30 minutes at 100x g;

5 8. Pipetted 43.61 μ l of Master Mix into PCR tubes;

9. Added appropriate reverse primer and DNA or water to
create the following primer pairs: (1, 2) + 5 μ l DNA; (1, 2) + 1 μ l; (1, 3) + 5 μ l DNA; (1, 3)
+ 1 μ l DNA; (1, 4) + 5 μ l DNA; (1, 4) + 1 μ l DNA; (1, 5) + 5 μ l DNA; (1, 5) + 1 μ l DNA;
(1, 6) + 5 μ l DNA; (1, 6) + 1 μ l DNA; (5, 8) + 5 μ l DNA; (7, 8) + 1 μ l DNA; (7, 9) + 5 μ l
20 DNA; (7, 9) + 1 μ l DNA; and (1, 2) + 5 (1, 4) + 5 μ l water;

10. Ran PCR;

11. Ran 1 % Agarose gel at 100 V for about 20 min.

Results: All E. coli primers showed amplification of target sequences, regardless of size.

Example 18

5 Purpose: To investigate the effects of 50 kGy irradiation on samples containing E. coli.

Materials:

1. E. coli spiked tendon (irradiated to 0 kGy and 50 kGy) +
6.65x10¹⁰ CFU/μl;
2. Taq PCR Core Kit (Qiagen, cat. no. 201223);
3. ProofStart Taq Polymerase (Qiagen, cat. no. 202205);
4. Primers: Ribo 7 and 8, and Ribo13, 14 and 15 SEQ ID
NOS.: 72, 73, and 74, respectively);
5. Agarose;
6. TAE Buffer;
7. Ethidium Bromide;
8. Microcon Centrifugal Filter Unit.

15 Procedure: 1. Prepared Master Mix as follows:

| | |
|--------------|----------|
| 10x Buffer | 60 μl |
| dNTP | 18 μl |
| pA (forward) | 16.68 μl |
| Taq | 12 μl |

ProofStart 12 µl

Water 452.64 µl;

2. Placed in Microcon and centrifuged for 30 minutes at 100x g;
3. Pipetted 47-61 µl master mix into each of 9 PCR tubes;
- 5 4. Added 1.39 µl of reverse primer and 1 µl DNA into appropriate tubes;
5. Ran PCR.
6. Ran 1% Agarose gel (8 µl sample + 1 µl sample buffer) at 100 V for about 20 minutes.

0 Results: Samples irradiated to 50 kGy showed progressive disappearance of bands with increasing amplicon size, indicating degradation of the E. coli genetic material caused by irradiation.

Example 19

15 Purpose: To show functionality of Mt-DNA primers for RT-PCR using large target amplicons.

Materials:

1. Tendon DNA (irradiated to 0 kGy and 50 kGy);
2. ROX 6 (1/10 dilution) molecular probes;
3. Primers: MITO 1, 2, 3, 4, and 5 (Figure 8) (SEQ ID NOS.: 90, 92, 95, 96, and 97, respectively);
- 20 4. MITO Probe 1 (Figure 8) (SEQ ID NO.: 91);

5. Human DNA;
6. Qiagen PCR Reagents;
7. Qiagen ProofStart.

Procedure:

1. Prepared the following mixtures:

| | | |
|---|----------------|---------------------------------|
| 5 | Buffer | 1.5 μ l |
| | dNTPs | 1.5 μ l |
| | MITO 1 | 2.5 μ l |
| | reverse primer | 2.5 μ l (MITO 2, 3, 4 or 5) |
| | MITO Probe | 2.5 μ l |
| 0 | Taq | 1 μ l |
| | PS | 1 μ l |
| | 1/10 ROX | 1 μ l |
| | water | 28 μ l |
| | DNA | 5 μ l |

2. Ran 40 PCR;
3. Ran 1% agarose gel (8 μ l product + 1 μ l sample loading buffer) at 100 V for about one hour.

Results: Mt-DNA primers were functional, regardless of amplicon size.

Example 20

20

Purpose: Real-time PCR amplification of human DNA (large amplicons).

Materials:

1. 10 ng of human control DNA; Calbiochem, Human Genomic DNA, Cat #HCD01, Lot # D10498;

2. Taq PCR Core Kit;

3. ProofStart DNA Polymerase;

5 4. Primers and Probes;

5. Agarose;

6. TAE;

7. Ethidium Bromide.

Procedure:

1. Prepared PCR Master Mix as follows:

| | | |
|----|-------------------------------------|-------------|
| 10 | 10x Buffer | 5 μ l |
| | dNTPs | 1.5 μ l |
| | Mito 1 (SEQ ID NO. 90) | 2.5 μ l |
| | Reverse Primer (Mito 5 or 7) | |
| | (SEQ ID NO. 97 or 99, respectively) | 2.5 μ l |
| 15 | MitoProbe 1 (SEQ ID NO. 91) | 2.5 μ l |
| | taq | 1 μ l |
| | Proof Start | 1 μ l |
| | water | 31 μ l |
| | water or DNA | 3 μ l |

20 2. Ran PCR (50 cycles);

3. Ran 8 µl PCR Products on 1% agarose gel and ran at 100 V
for about 20 minutes.

Results: Target sequences greater then 8,000 nucleic acid residues can be
5 successfully amplified with Real-time PCR.

Example 21

Purpose: Real-time PCR on fibular bone rings to detect bacterial contamination
in unirradiated bone samples.

0

Materials:

1. Bacterial extracts from bones;
2. Taq PCR Core Kit (Qiagen, Cat#201223);
3. ProofStart DNA Polymerase (Qiagen, Cat#202205);
4. Primers: Ribo 7 and 10 (SEQ ID NOS. 100 and 69,
.5 respectively);
5. Probe: FAM-RIBO (SEQ ID NO. 101);
6. Optically clear plates and seals;

Procedure: 1. Prepared PCR setup as follows:

| | | Per run | x23 |
|----|------------|---------|---------|
| 20 | 10x Buffer | 5 µl | 115 µl |
| | dNTPs | 1.5 µl | 34.5 µl |
| | pA | 3.5 µl | 80.5 µl |

| | | | |
|---|-------------|---------------|----------------|
| | pB | 3.5 μ l | 80.5 μ l |
| | Probe | 2.5 μ l | 57.5 μ l |
| | taq | 0.25 μ l | 5.75 μ l |
| | Proof Start | 1 μ l | 23 μ l |
| 5 | water | 30.75 μ l | 707.25 μ l |

2. Aliquot 48 μ l into A4-7, B4-7, C4-7, D4-7, and H11-12;
3. Pipet 2 μ l of appropriate DNA into each well
4. Seal plate and run "long" program on the thermocycler (40 cycles).

10 Results: Of 103 bone samples, 40% were found to be contaminated with bacteria.

Having now fully described this invention, it will be understood to those of ordinary skill in the art that the methods of the present invention can be carried out with a wide and equivalent range of conditions, formulations and other parameters without departing from the scope of the invention or any embodiments thereof.

15 All patents and publications cited herein are hereby fully incorporated by reference in their entirety. The citation of any publication is for its disclosure prior to the filing date and should not be construed as an admission that such publication is prior art or that the present invention is not entitled to antedate such publication by virtue of prior invention.

20 The foregoing embodiments and advantages are merely exemplary and are not to be construed as limiting the present invention. The present teaching can be readily applied to other types of apparatuses. The description of the present invention is intended to be

illustrative, and not to limit the scope of the claims. Many alternatives, modifications, and variations will be apparent to those skilled in the art. In the claims, means-plus-function clauses are intended to cover the structures described herein as performing the recited function and not only structural equivalents but also equivalent structures.

WHAT IS CLAIMED IS:

1. A method for analyzing a target nucleic acid sequence in a biological material, said method comprising:

5 (i) adding to said biological material an effective amount of at least two nucleic acid primers,

wherein said nucleic acid primers hybridize under stringent conditions to predetermined nucleic acid sequences of said target nucleic acid sequence that are separated by at least about 750 nucleic acid residues,

10 (ii) amplifying said target nucleic acid sequence by polymerase chain reaction, said polymerase chain reaction comprising adding a polymerase to said biological material and primers to form an amplification mixture and thermally cycling said amplification mixture between at least one denaturation temperature and at least one elongation temperature,

15 wherein said elongation temperature is not more than about 70°C and said denaturation temperature is not more than about 95°C, and further wherein during each thermal cycle said amplification mixture is maintained at said denaturation temperature for a period of not more than about 30 seconds and at said elongation temperature for a period of not less than about 1 minute; and

20 (iii) detecting and quantifying said target nucleic acid sequence.

2. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by at least about 1000 nucleic acid residues of said target nucleic acid sequence

3. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by at least about 2000 nucleic acid residues of said target nucleic acid sequence

5

4. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by at least about 3000 nucleic acid residues of said target nucleic acid sequence

.0

5. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by at least about 4000 nucleic acid residues of said target nucleic acid sequence.

.5

6. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by at least about 5000 nucleic acid residues of said target nucleic acid sequence.

20

7. The method according to claim 1, wherein said predetermined nucleic acid sequences of said target nucleic acid sequence are separated by only about 500 nucleic acid residues of said target nucleic acid sequence.

8. The method according to claim 1, wherein said step (i) further comprises adding at least one nucleic acid probe to said biological material.

25

9. The method according to claim 8, wherein said nucleic acid probe is selected from the group consisting of 5' nuclease probes, hairpin probes, adjacent probes, sunrise probes and scorpion probes.

30

10. The method according to claim 1, wherein said elongation temperature is between about 60°C and about 69°C.

11. The method according to claim 1, wherein said elongation temperature is between about 65°C and about 69°C.

5 12. The method according to claim 1, wherein said denaturation temperature is between about 90°C and about 95°C.

13. The method according to claim 1, wherein said denaturation temperature is between about 93°C and about 95°C.

0

14. The method according to claim 1, wherein during each thermal cycle said amplification mixture is maintained at said denaturation temperature for a period of not more than about 20 seconds.

5 15. The method according to claim 1, wherein during each thermal cycle said amplification mixture is maintained at said denaturation temperature for a period of not more than about 10 seconds.

16. The method according to claim 1, wherein during each thermal cycle said
0 amplification mixture is maintained at said elongation temperature for a period of not less than about 2 minutes.

17. The method according to claim 1, wherein during each thermal cycle said
5 amplification mixture is maintained at said elongation temperature for a period of not less than about 3 minutes.

18. The method according to claim 1, wherein the period during which said
amplification mixture is maintained at said elongation temperature during each thermal cycle
is increased by a period of about 5 seconds for each successive thermal cycle.

0

19. The method according to claim 1, wherein said amplification mixture is thermally cycled for at least 30 cycles.

20. The method according to claim 1, wherein said amplification mixture is
5 thermally cycled for at least 40 cycles.

21. The method according to claim 1, wherein said amplification mixture is thermally cycled for at least 50 cycles.

22. The method according to claim 1, wherein said biological material has been
10 subjected to an environment or process that may have altered said target nucleic acid sequence.

23. The method according to claim 1, wherein said polymerase is a *Taq*
15 polymerase.

24. The method according to claim 1, wherein said polymerase is a proof-reading
Taq polymerase.

25. The method according to claim 1, wherein said amplification mixture further
20 comprises at least one thermostable inorganic pyrophosphatase.

26. The method according to claim 25, wherein the ratio of *Taq* polymerase to
thermostable inorganic pyrophosphatase is about 5:1.
25

CI0042PCTseqlisting.ST25
SEQUENCE LISTING

<110> Clearant, Inc.
McKenney, Keith
Gillmeister, Lidja
Marlowe, Kristina
Armistad, David

<120> Real-Time Polymerase Chain Reaction Using Large Target Amplicons

<130> CI-0042PCT

<150> 10/361,004
<151> 2003-02-10

<160> 101

<170> PatentIn version 3.2

<210> 1
<211> 5594
<212> DNA
<213> B19 virus

<400> 1
 ccaaatcaga tgccgccggt cgccgccggt aggcgggact tccggtacaa gatggcggac 60
 aattacgtca tttcctgtga cgtcatttcc tgtgacgtca cttccggtgg gcgggacttc 120
 cggaattagg gttggctctg ggccagcttg cttggggttg cttgacact aagacaagcg 180
 gcgcgccgct tgtcttagtg gcacgtcaac cccaagcgct ggcccagagc caaccctaata 240
 tccggaagtc ccgccaccg gaagtgcagt cacaggaaat gacgtcacag gaaatgacgt 300
 aattgtccgc catcttgtac cggaagtccc gcctaccggc ggcgaccggc ggcatctgat 360
 ttggtgtctt cttttaaatt ttagcgggct tttttccgc cttatgcaaa tgggcagcca 420
 ttttaagtgt ttactataa ttttattggt cagttttgta acggttaaaa tgggcggagc 480
 gtaggcgggg actacagtat atatagcacg gactgcccgc agctctttct ttctgggctg 540
 ctttttcctg gactttcttg ctgttttttg tgagctaact aacagggtatt tatactactt 600
 gttaacatac taacatggag ctatttagag ggtgtcttca agtttcttct aatgttctgg 660
 actgtgctaa cgataactgg tgggtgctctt tactggattt agacacttct gactgggaac 720
 cactaactca tactaacaga ctaatggcaa tatacttaag cagtgtggct tctaagcttg 780
 actttaccgg ggggccacta gcggggtgct tgtacttttt tcaagtagaa tgtaacaaat 840
 ttgaagaagg ctatcatatt catgtggtta ttggggggcc aggggttaaac ccagaaacc 900
 tcacagtgtg tgtagagggg ttatttaata atgtacttta tcaccttgta actgaaaatg 960
 taaagctaaa atttttgcc ggaatgacta caaaaggcaa atactttaga gatggagagc 1020
 agtttataga aaactattta atgaaaaaaa tacctttaaa tgttgatgg tgtgttacta 1080
 atattgatgg atatataga' acctgtat' ctgctacttt tagaagggga gcttgccatg 1140
 ccaagaaacc ccgcattacc acagccataa atgacactag tagtgatgct ggggagtcta 1200
 gcggcacagg ggcagaggtt gtgccaat' atgggaaggg aactaaggct agcataaagt 1260
 ttcaaactat ggtaaactgg ttgtgtgaaa acagagtgtt tacagaggat aagtggaaac 1320

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|------------|-------------|-------------|-------------|------------|------|
| tagttgactt | taaccagtac | actttactaa | gcagtagtca | cagtgggaagt | tttcaaattc | 1380 |
| aaagtgcact | aaaactagca | atttataaag | caactaattt | agtgcctaca | agcacatttc | 1440 |
| tattgcatac | agactttgag | caggttatgt | gtattaaaga | caataaaatt | gttaaattgt | 1500 |
| tactttgtca | aaactatgac | cccctattag | tggggcagca | tgtgttaaag | tggattgata | 1560 |
| aaaaatgtgg | caagaaaaat | acactgtggt | tttatgggcc | gccaagtaca | ggaaaaacaa | 1620 |
| acttggcaat | ggccattgct | aaaagtgttc | cagtatatgg | catggttaac | tggaataatg | 1680 |
| aaaactttcc | atttaatgat | gtagcagggg | aaagcttggt | ggtctgggat | gaaggattaa | 1740 |
| ttaagtctac | aattgtagaa | gctgcaaaag | ccattttagg | cgggcaaccc | accagggtag | 1800 |
| atcaaaaaat | gcgtggaagt | gtagctgtgc | ctggagtacc | tgtggttata | accagcaatg | 1860 |
| gtgacattac | ttttgttgta | agcgggaaca | ctacaacaac | tgtacatgct | aaagccttaa | 1920 |
| aagagcgaat | ggtaaagtta | aactttactg | taagatgcag | ccctgacatg | gggttactaa | 1980 |
| cagaggctga | tgtacaacag | tggcttacat | ggtgtaatgc | acaaagctgg | gaccactatg | 2040 |
| aaaactgggc | aataaactac | acttttgatt | tccctggaat | taatgcagat | gccctccacc | 2100 |
| cagacctcca | aaccacccca | attgtcacag | acaccagtat | cagcagcagt | ggtggtgaaa | 2160 |
| gctctgaaga | actcagtga | agcagctttt | ttaacctcat | caccccaggc | gcctggaaca | 2220 |
| ctgaaacccc | gcgctctagt | acgcccattc | ccgggaccag | ttcaggagaa | tcatttgtcg | 2280 |
| gaagctcagt | ttcctccgaa | gttgtagctg | catcgtggga | agaagccttc | tacacacctt | 2340 |
| tggcagacca | gtttcgtgaa | ctgttagttg | gggttgatta | tgtgtgggac | ggtgtaaggg | 2400 |
| gtttacctgt | gtgttggtg | caacatatta | acaatagtgg | gggaggcttg | ggactttgtc | 2460 |
| cccattgcat | taatgtaggg | gcttgggtata | atggatggaa | atttcgagaa | tttaccacag | 2520 |
| atttggtgcg | gtgtagctgc | catgtgggag | cttctaattc | cttttctgtg | ctaacctgca | 2580 |
| aaaaatgtgc | ttacctgtct | ggattgcaaa | gctttgtaga | ttatgagtaa | agaaagtggc | 2640 |
| aaatgggtggg | aaagtgatga | taaatttgct | aaagctgtgt | atcagcaatt | tgtggaattt | 2700 |
| tatgaaaagg | ttactggaac | agacttagag | cttattcaaa | tattaaaaga | tcactataat | 2760 |
| atttcttttag | ataatcccct | agaaaaccca | tcctctctgt | ttgacttagt | tgctcgtatt | 2820 |
| aaaaataacc | ttaaaaactc | tccagactta | tatagtcatc | attttcaaag | tcatggacag | 2880 |
| ttatctgacc | acccccatgc | cttatcatcc | agtagcagtc | atgcagaacc | tagaggagaa | 2940 |
| aatgcagtat | tatctagtga | agacttacac | aagcctgggc | aagttagcgt | acaactaccc | 3000 |
| ggtactaact | atgttggggc | tggcaatgag | ctacaagctg | ggcccccgca | aagtgtgttt | 3060 |
| gacagtgctg | caaggattca | tgactttagg | tatagccaac | tggctaagtt | gggaataaat | 3120 |
| ccatatactc | attggactgt | agcagatgaa | gagcttttaa | aaaatataaa | aaatgaaact | 3180 |
| gggtttcaag | cacaagtagt | aaaagactac | tttactttaa | aaggtgcagc | tgcccctgtg | 3240 |
| gcccattttc | aaggaagttt | gccggaagtt | cccgtttaca | acgcctcaga | aaaataccca | 3300 |
| agcatgactt | cagttaattc | tgcaagaagc | agcactgggtg | caggaggggg | tggcagtaat | 3360 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|------------|------------|-------------|------------|------------|-------------|------|
| cctgtcaaaa | gcatgtggag | tgagggggcc | acttttagtg | ccaactctgt | aacttgtaca | 3420 |
| ttttccagac | agtttttaat | tccttatgac | ccagagcacc | attataaggt | gttttctccc | 3480 |
| gcagcaagca | gctgccacaa | tgccagtgga | aaggaggcaa | aggtttgac | aattagtccc | 3540 |
| ataatgggat | actcaacccc | atggagatat | ttagatttta | atgctttaaa | tttatttttt | 3600 |
| tcacctttag | agtttcagca | cttaattgaa | aattatggaa | gtatagctcc | tgatgcttta | 3660 |
| actgtaacca | tatcagaaat | tgctgttaag | gatgttacag | acaaaactgg | agggggggta | 3720 |
| caggttactg | acagcactac | agggcgcccta | tccatgttag | tagaccatga | atacaagtac | 3780 |
| ccatatgtgt | taggacaagg | tcaggatact | ttagccccag | aacttcctat | ttgggtatac | 3840 |
| tttccccctc | aatatgctta | cttaacagta | ggagatgtta | acacacaagg | aatctctgga | 3900 |
| gacagcaaaa | aattagcaag | tgaagaatca | gcattttatg | ttttggaaca | cagttctttt | 3960 |
| cagcttttag | gtacaggagg | tacagcaact | atgtcttata | agtttctctc | agtgtcccca | 4020 |
| gaaaatttag | agggtcgcag | tcaacacttt | tatgaaatgt | acaatccctt | atacggatcc | 4080 |
| cgcttagggg | ttcctgacac | attaggagggt | gacccaaaat | ttagatcttt | aacacatgaa | 4140 |
| gaccatgcaa | ttcagcccca | aaacttcacg | ccagggccac | tagtaaaact | agtgtctaca | 4200 |
| aaggagggag | acagctctaa | tactggagct | ggaaaagcct | taacaggcct | tagcacaggc | 4260 |
| acctctcaaa | acactagaat | atccttacgc | cctgggccag | tgtcacagcc | ataccaccac | 4320 |
| tgggacacag | ataaatatgt | tccaggaata | aatgccattt | ctcatggtca | gaccacttat | 4380 |
| ggtaacgctg | aagacaaaga | gtatcagcaa | ggagtgggta | gattttcaaa | tgaaaaagaa | 4440 |
| cagctaaaac | agttacaggg | tttaaacatg | cacacctatt | tccccaataa | aggaaccag | 4500 |
| caatatacag | atcaaattga | gcgcccccta | atggtgggtt | ctgtatggaa | cagaagagcc | 4560 |
| cttactatg | aaagccagct | gtggagtaaa | attccaaatt | tagatgacag | ttttaaaact | 4620 |
| cagtttgacg | ccttaggagg | atggggtttg | catcagccac | ctcctcaa | atttttaaaa | 4680 |
| atattaccac | aaagtgggcc | aattggagggt | attaaatcaa | tgggaattac | tacctagtt | 4740 |
| cagtatgccg | tgggaattat | gacagtaact | atgacattta | aattggggcc | ccgtaaagct | 4800 |
| acgggacggt | ggaatcctca | acctggagta | tatccccgc | acgcagcagg | tcattttacca | 4860 |
| tatgtactat | atgacccac | agctacagat | gcaaaacaac | accacaggca | tggatacgaa | 4920 |
| aagcctgaag | aattgtggac | agccaaaagc | cgtgtgcacc | cattgtaaac | actccccacc | 4980 |
| gtgccctcag | ccaggatgcg | taactaaacg | cccaccagta | ccaccagac | tgtacctgcc | 5040 |
| ccctcctgta | cctataagac | agcctaacac | aaaagatata | gacaatgtag | aatttaagta | 5100 |
| cttaaccaga | tatgaacaac | atgttattag | aatgttaaga | ttgtgtaata | tgtatcaaaa | 5160 |
| tttagaaaaa | taaacatttg | ttgtgggtta | aaaattatgt | tggtgcgctt | taaaaattta | 5220 |
| aaagaagaca | ccaaatcaga | tgccgccggt | cgccgccggt | aggcgggact | tccggtacaa | 5280 |
| gatggcggac | aattacgtca | tttccctgtga | cgtcatttcc | tgtgacgtca | cttccggtgg | 5340 |
| gcgggacttc | cggaattagg | gttggtctctg | ggccagcgct | tgggggtgac | gtgccactaa | 5400 |

CI0042PCTseqlisting.ST25

gacaagcggc gcgccgcttg tcttagtgtc aaggcaaccc caagcaagct ggcccagagc 5460
 caaccctaatt tccggaagtc ccgcccaccg gaagtgcgacg cacaggaaat gacgtcacag 5520
 gaaatgacgt aattgtccgc catcttgtag cggaagtccc gcctaccggc ggcgaccggc 5580
 ggcatctgat ttgg 5594

<210> 2
 <211> 3221
 <212> DNA
 <213> Hepatitis B virus

<400> 2
 ttccactgcc ttccaccaag ctctgcaaga ccccagagtc aggggtctgt attttcctgc 60
 tgggtggctcc agttcaggaa cagtaaaccg tgctccgaat attgcctctc acatctcgtc 120
 aatctccgcg aggaccgggg accctgtgac gaacatggag aacatcacat caggattcct 180
 aggacccctg cccgtgttac aggcgggggt tttcttggtg acaagaatcc tcacaatacc 240
 gcagagtcta gactcgtggg ggacttctct caattttcta gggggatcac ccgtgtgtct 300
 tggccaaaat tcgcgatccc caacctccaa tcactacca acctcctgtc ctccaatttg 360
 tcctggttat cgctggatgt gtctgcggcg ttttatcata ttcctcttca tcctgctgct 420
 atgcctcatc ttcttattgg ttcttctgga ttatcaagg atgttgccc tttgtcctct 480
 aattctagga tcaacaacaa ccagtacggg accatgcaaa acctgcacga ctctgctca 540
 aggcaactct atgtttccct catgttgctg tacaaaacct acggatggaa attgcacctg 600
 tattcccatc ccatcgtctt gggctttcgc aaaataccta tgggagtggg cctcagtcg 660
 tttctcttgg ctcagtttac tagtgccatt tggtcagtgg ttcgtagggc tttccccac 720
 tgtttggtct tcagctatat ggatgatgtg gtattggggg ccaagtctgt acagcatcgt 780
 gagttccttt ataccgctgt taccaatttt cttttgtctc tgggtataca tttaaacctt 840
 aacaaaacaa aaagatgggg ttattcccta aacttcatgg gttatgtaat tgggaagtgg 900
 ggaacattgc cacaggatca tattgtacaa aaaatcaaac actgttttag aaaacttctt 960
 gttaacaggc ctattgattg gaaagtatgt caaagaattg tgggtctttt gggctttgct 1020
 gtccttttta cacaatgtgg atatcctgcc ttaatgccct tgtatgcatg tatacaagct 1080
 aaacaggctt tcactttctc gccaaacttac aaggcctttc taagtaaaca gtacatgaac 1140
 ctttaccctg ttgtcggca acggcctggt ctgtgccaa gatttgctga tgcaaccccc 1200
 actggctggg gcttggccat aggcacatcag cgcatgcgcg gaacctttgt ggctcctctg 1260
 ccgatccata ctgcggaact cctagccgct tgttttgctc gcagccgggtc tggagcgaaa 1320
 ctcatcgga ctgacaattc tgcgtcctc tcgcggaaat atacctcgtt tccatggcta 1380
 ctaggctgtg ctgccaaact gatccttcgc gggacgtcct ttgtttacgt cccgtcggcg 1440
 ctgaatcccg cggacgacct ctctcggggc cgcttgggac tctctcgtcc ccttctccgt 1500
 ctgccgttcc agccgaccac ggggcgcacc tctctttacg cggctctccc gtctgtgcct 1560

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|-------------|------------|------------|------------|------------|------|
| tctcatctgc | cgggccgtgt | gcacttcgct | tcacctctgc | acgttgcatg | gagaccaccg | 1620 |
| tgaacgcca | tcagatcctg | ccaaggtct | tacataagag | gactcttgga | ctcccagcaa | 1680 |
| tgtcaacgac | cgaccttgag | gcctacttca | aagactgtgt | gtttaaggac | tgggaggagc | 1740 |
| tgggggagga | gattaggtta | aaggtctttg | tattaggagg | ctgtaggcac | aaattggtct | 1800 |
| gcgcaccagc | accatgcaac | tttttcacct | ctgcctaata | atctcttgta | catgtccac | 1860 |
| tgttcaagcc | tccaagctgt | gccttgggtg | gctttggggc | atggacattg | acccttataa | 1920 |
| agaatttgga | gctactgtgg | agttactctc | gtttttgcct | tctgacttct | ttccttcctg | 1980 |
| cagagatctc | ctagacaccg | cctcggtctc | gtatcgggaa | gccttagagt | ctcctgagca | 2040 |
| ttgctcacct | caccataaccg | cactcaggca | agccattctc | tgctgggggg | aattgatgac | 2100 |
| tctagctacc | tgggtgggta | ataatttgga | agatccagca | tccagggatc | tagtagtcaa | 2160 |
| ttatgttaat | actaacatgg | gattaaagat | caggcaactc | ttgtggtttc | atatctcttg | 2220 |
| ccttactttt | ggaagagaaa | ctgtacttga | atatttggtc | tctttcggag | tgtggattcg | 2280 |
| cactcctcca | gcctatagac | caccaaatac | ccctatctta | tcaacacttc | cggaaactac | 2340 |
| tggtgttaga | cgacgggacc | gaggcaggtc | ccctagaaga | agaactccct | cgcctcgag | 2400 |
| acgcagatct | caatcgccgc | gtcgcagaag | atctcaatct | cgggaatctc | aatgttagta | 2460 |
| ttccttgagc | tcataagggtg | ggaaacttca | ctgggcttta | ttcctctaca | gcacctatct | 2520 |
| ttaatcctga | atggcaaact | ccttcctttc | ctaaaattca | tttacaagag | gacattatta | 2580 |
| atagggtgtca | acaatttggtg | ggccctctca | ctgtaaatga | aaagagaaga | ttgaaattaa | 2640 |
| ttatgcctgc | tagattctat | cctaccacac | ctaaatattt | gcccttagac | aaaggaatta | 2700 |
| aaccttatta | tccagatcag | gtagttaatc | attacttcca | aaccagacat | tatttacata | 2760 |
| ctcttttgaa | ggcgggtatt | ctatataaga | gagaaaccac | acgtagcgca | tcattttgag | 2820 |
| ggtcaccata | ttcttgggaa | caagagctac | agcatgggag | gttggtcatc | aaaacctcgc | 2880 |
| aaaggcatgg | ggacgaatct | ttctgttccc | aacctcttgg | gattctttcc | cgatcatcag | 2940 |
| ttggaccctg | tattcggagc | caactcaaac | aatccagatt | gggacttcaa | cccatcaag | 3000 |
| gaccactggc | cagcagccaa | ccaggtagga | gtgggagcat | tcgggccagg | gttcacccct | 3060 |
| ccacacggcg | gtgttttggtg | gtggagccct | caggctcagg | gcatgttgac | cccagtgta | 3120 |
| acaattcctc | ctcctgcctc | cgccaatcgg | cagtcaggaa | ggcagcctac | tcccatctct | 3180 |
| ccacctctaa | gagacagtca | tcctcaggcc | atgcagtgga | a | | 3221 |

<210> 3
 <211> 5075
 <212> DNA
 <213> Porcine parvovirus

| | | |
|------------|------------|------------|
| <400> 3 | | |
| aatcttttaa | ctgaccaact | gtctttgcgt |
| atggtgacgt | gatgacgcgc | gctacgcgcg |
| | | 60 |
| ctgccttcgg | cagtcacacg | tcaccatcag |
| caaagacagt | tggtcagttt | aaagattaat |
| | | 120 |
| aagacattcc | attggctgaa | aagaggcggg |
| aaattcaaaa | aaagaggcgg | gaaaaaaaga |
| | | 180 |

CI0042PCTseqlisting.ST25

| | |
|---|------|
| ggtggagcct aacactataa atacagttgc ttacttcagt tagttccttt ctgcttcaga | 240 |
| ctgcacttcg ctccagagac acagctacaa actactctca gctactgcag catggcagcg | 300 |
| ggaaacactt actcggaga ggtactaaaa gctaccaact ggcttcaaga taatgctcaa | 360 |
| aaagaagcat tctcttatgt atttaaaaca caaaaagtca atctaaatgg aaaagaaatt | 420 |
| gcttgggaata actacaacaa agatacaaca gatgcggaaa tgataaacct acaaagagga | 480 |
| gcagaaacat catgggacca ggcaacagac atggaatggg aatcagaaat cgacagcctc | 540 |
| acaaaacggc aagtactgat ttttgactct cttgttaaaa aatgtctctt tgaaggata | 600 |
| ttgcaaaaga acctaagtcc aagtgactgc tactggttca tacagcatga acatggtcaa | 660 |
| gatactggct atcactgcc tgtactacta ggtggaaaag gcttacaaca agcaatggga | 720 |
| aaatgggttca gaaaacaatt aaacaattta tggagtagat ggttaataat gcaatgcaaa | 780 |
| gtacctctaa caccagttga aagaataaaa ttaagggaaat tagcagagga tggtagtggtg | 840 |
| gtatcgctac taacctacac tcacaaacaa actaaaaaac aatatacaaa aatgactcat | 900 |
| tttggaataa tgattgctta ctacttccta aataaaaaaa gaaagacaac tgaaagagag | 960 |
| catggatatt atctcagctc agattctggc ttcattgacaa atttcttaaa agaaggcgag | 1020 |
| agacacttag tcagtcacct atttactgaa gcaataaac ctgaaactgt ggaaacaacg | 1080 |
| gttactacag ctcaaggaagc caaaagaggc agaatacaaa caaaaaaaga agtaagcata | 1140 |
| aaatgcacaa taagagactt ggttaataaa agatgtacta gcatagaaga ctggatgatg | 1200 |
| acagatccag acagttatat agaaatgatg gctcaaaccg gaggagaaaa tttaatcaaa | 1260 |
| aatacactag aaataacaac tcttactcta gcaagaacaa aaacagcata tgacttaata | 1320 |
| cttgaaaagg caaaaccaag catgctacca acatttaata ttagcaatac aagaacatgt | 1380 |
| aaaatattca gcatgcacaa ttggaactac attaaagtct gccatgctat aacttggtga | 1440 |
| ctaaacagac aaggaggaaa aagaaataca attctatttc atgggccagc atcaacagga | 1500 |
| aaaagtataa ttgctcaaca cattgcaaac ttagttggta atggttggtg ctacaatgca | 1560 |
| gccaatgtga actttccatt taatgactgt acaataaaaa acttaatatg gattgaagaa | 1620 |
| gcaggaaact tctctaacca agtaaacc aa ttcaaagcca tatgttcagg tcaaacaatt | 1680 |
| agaattgacc aaaaaggtaa aggaagcaaa caaattgaac caactcctgt aataatgact | 1740 |
| acaaatgaag acataactaa agttagaata ggatgagagg aaagaccaga acatacacia | 1800 |
| ccaataagag acagaatggt aaacataaac ctaaccagaa aactgccagg tgattttgga | 1860 |
| cttttagaag aaactgaatg gccactaata tgtgcttggt tggtaaagaa aggttacaa | 1920 |
| gcaacaatgg ctagctatat gcatcattgg ggaaatgtac ctgattggtc agaaaaatgg | 1980 |
| gaggagccaa aaatgcaaac cccaataaat acaccaacag actctcagat ttccacatca | 2040 |
| gtgaaaactt cgccagcgga caacaactac gcagcaactc caatacagga ggacctggat | 2100 |
| ttagcttttag ccttgagacc gtggagcgag ccaacaacac caactttcac caactgcac | 2160 |
| ttaactccaa caccgccaga ttcagcaata cggacaccaa gtccaacttg gtcggaaata | 2220 |

CI0042PCTseqlisting.ST25

| | |
|--|------|
| gaaaccgaca taagagcctg ctttggtgaa aactgtgcac ccacaacaaa ccttgaataa | 2280 |
| ggtaggatgg cgcctcctgc aaaaagagca agaggtaagg gtagttttta gggggtggtg | 2340 |
| ggcatacata taaaactaac tgcaaataat ttttttatat attacaggac taactctacc | 2400 |
| aggatacaaa taccttggtc caggaaactc actagaccaa ggagaaccaa ctaatccatc | 2460 |
| agacgccgca gcaaaagaac acgacgaagc ctacgacaaa tacataaaat ctggaaaaaa | 2520 |
| tccatacttc tactttctcag cagctgatga aaaattcata aaagaaactg aacacgcaaa | 2580 |
| agactacgga ggtaaaattg gacattactt cttcagagca aagcgtgcct ttgctccaaa | 2640 |
| actctcagaa acagactcac caactacatc tcaacaacca gaggtaagaa gatcgccgag | 2700 |
| aaaacaccca ggggtctaac caccaggaag aagacctgct ccaagacata tttttataaa | 2760 |
| cttagctaaa aaaaaagcta aagggaacatc taatacaaac tctaactcaa tgagtgaaaa | 2820 |
| tgtggaacaa cacaacccta ttaatgcagg cactgaattg tctgcaacag gaaatgaatc | 2880 |
| tgggggtggg ggcggcggtg gcgggggtag ggggtgctggg ggggttggtg tgtctacagg | 2940 |
| tactttcaat aatcaaacag aattttcaata cttgggggag ggcttggtta gaatcactgc | 3000 |
| acacgcatca agactcatac atctaaatat gccagaacac gaaacataca aaagaatata | 3060 |
| tgtactaaat tcagaatcag ggggtggcggg acaaatggta caagacgatg cacacacaca | 3120 |
| aatggttaaca ccttggtcac taatagatgc taacgcatgg ggagtgtggt tcaatccagc | 3180 |
| ggactggcag ttaatatcca acaacatgac agaaataaac ttagttagtt ttgaacaaga | 3240 |
| aatattcaat gtagtactta aaacaattac agaatcagca acctcaccac caacaaaat | 3300 |
| atataataat gatctaactg caagcttaat ggctgcacta gacaccaata acacacttcc | 3360 |
| atacacacca gcagcaccta gaagtgaac acttggtttt tatccatggt tacctacaaa | 3420 |
| accaactcaa tacagatatt acctatcatg catcagaaac ctaaatccac caacatacac | 3480 |
| tggacaatca caacaaataa cagactcaat acaaacagga ctacacagtg acattatggt | 3540 |
| ctacacaata gaaatgcag taccaattca tcttctaaga acaggagatg aattctccac | 3600 |
| aggaatatat cactttgaca caaaaccact aaaattaact cactcatggc aaacaaacag | 3660 |
| atctctagga ctgcctccaa aactactaac tgaacctacc acagaaggag accaacaccc | 3720 |
| aggaacacta ccagcagcta acacaagaaa aggttatcac caaacaatta ataatagcta | 3780 |
| cacagaagca acagcaatta ggccagctca ggtaggatat aatacaccat acatgaattt | 3840 |
| tgaatactcc aatggtggac cattttctaac tcctatagta ccaacagcag acacacaata | 3900 |
| taatgatgat gaaccaaagt gtgctataag atttacaatg gattaccaac atggacactt | 3960 |
| aaccacatct tcacaagagc tagaaagata cacattcaat ccacaaagta aatgtggaag | 4020 |
| agctccaaag caacaattta atcaacaggc accactaaac ctagaaaata caaataatgg | 4080 |
| aacactttta ccttcagatc caataggagg gaaatctaac atgcatttca tgaatacact | 4140 |
| caatacatat ggaccattaa cagcactaaa caatactgca cctgtatttc caaatgggtca | 4200 |
| aatatgggat aaagaacttg atacagatct aaaacctaga ctacatgtta cagctccatt | 4260 |

CI0042PCTseqlisting.ST25

tgtttgtaaa aacaatccac caggacaact atttgtaaaa atagcaccaa acctaacaga 4320
 tgatttcaat gctgactctc ctcaacaacc tagaataata acttattcaa acttttggtg 4380
 gaaaggaaca ctaacattca cagcaaaaat gagatccagt aatatgtgga accctattca 4440
 acaacacaca acaacagcag aaaacattgg taactatatt cctacaaata ttggtggcat 4500
 aagaatgttt ccagaatatt cacaacttat accaagaaaa ttataactaga aataactctg 4560
 taaataaaaa ctcagttact tggttaatca tgtactacta tcattgtata cttcaataaa 4620
 aataaattgt aaaatcaata aaactaagtt acttagtttc tgtataccta tactagaaat 4680
 aactctgtaa ataaaaactc agttacttgg ttaatcatgt actactatca ttgtatactt 4740
 caataaaaaat aaattgtaaa atcaataaaa ctaagttact tagtttctgt ataccaatta 4800
 tcccaaaaaa acaataaaaat tttaaaaaga aacaagctct catgtgttta ctattaacta 4860
 aaccaaccac acttatatga ccttatgtct ttaggggtggg tgggtgggaa ttactatgta 4920
 ttcttttgag ttagttgggc gcctttgggc gactaaccaa gcggctctgc cgcttggtta 4980
 gtcgcacggc gaccaactaa ctcaaaggaa tacatagtaa tccccacca cccaccctaa 5040
 agacataagg tcatataagt gtggttggtt tagtt 5075

<210> 4
 <211> 11703
 <212> DNA
 <213> Sindbis virus

<400> 4
 attgacggcg tagtacacac tattgaatca aacagccgac caattgcact accatcacaa 60
 tggagaagcc agtagtaaac gtagacgtag acccccagag tccgtttgtc gtgcaactgc 120
 aaaaaagctt cccgcaattt gaggtagtag cacagcaggt cactccaaat gaccatgcta 180
 atgccagagc attttcgcat ctggccagta aactaatcga gctggagggt cctaccacag 240
 cgacgatctt ggacataggc agcgcaccgg ctcttagaat gttttccgag caccagtatc 300
 attgtgtctg ccccatgcgt agtccagaag acccggaccg catgatgaaa tacgccagta 360
 aactggcgga aaaagcgtgc aagattacaa acaagaactt gcatgagaag attaaggatc 420
 tccggaccgt acttgatacg ccggatgctg aaacaccatc gctctgcttt cacaacgatg 480
 ttacctgcaa catgcgtgcc gaatattccg tcatgcagga cgtgtatatc aacgctcccg 540
 gaactatcta tcatcaggct atgaaaggcg tgcggaccct gtactggatt ggcttcgaca 600
 cccccagtt catgttctcg gctatggcag gttcgtaccc tgcgtacaac accaactggg 660
 ccgacgagaa agtccttgaa gcgcgtaaca tcggactttg cagcaciaag ctgagtgaag 720
 gtaggacagg aaaattgtcg ataattgagga agaaggagtt gaagcccggg tcgcggtttt 780
 atttctccgt aggatcgaca ctttatccag aacacagagc cagcttgag agctggcatc 840
 ttccatcggg gttccacttg aatggaaagc agtcgtacac ttgccgctgt gatacagtgg 900
 tgagttgcga aggctacgta gtgaagaaaa tcaccatcag tcccgggatc acgggagaaa 960

CI0042PCTseqlisting.ST25

| | |
|---|------|
| ccgtgggata cgcggttaca cacaatagcg agggcttctt gctatgcaaa gttactgaca | 1020 |
| cagtaaaagg agaacgggta tcgttccctg tgtgcacgta catcccggcc accatatgcg | 1080 |
| atcagatgac tggataatg gccacggata tatcacctga cgatgcacaa aaacttctgg | 1140 |
| ttgggctcaa ccagcgaatt gtcattaacg gtaggactaa caggaacacc aacaccatgc | 1200 |
| aaaattacct tctgccgac atagcacaa gggttcagca atgggctaag gagcgcaagg | 1260 |
| atgatcttga taacgagaaa atgctgggta ctagagaacg caagcttacg tatggctgct | 1320 |
| tgtgggcgtt tcgcactaag aaagtacatt cgttttatcg cccacctgga acgcagacct | 1380 |
| gcgtaaaagt cccagcctct tttagcgctt ttcccatgtc gtccgtatgg acgacctctt | 1440 |
| tgcccatgtc gctgaggcag aaattgaaac tggcattgca accaaagaag gagggaaaaac | 1500 |
| tgctgcaggt ctcgaggagaa ttagtcatgg aggccaaaggc tgcttttgag gatgctcagg | 1560 |
| aggaagccag agcggagaag ctccgagaag cacttccacc attagtggca gacaaaggca | 1620 |
| tcgaggcagc cgcagaagtt gtctgcgaag tggaggggct ccaggcggaac atcggagcag | 1680 |
| cattagttga aaccccgcg ggtcacgtaa ggataatacc tcaagcaaat gaccgtatga | 1740 |
| tcggacagta tatcgttgtc tcgccaaact ctgtgctgaa gaatgccaaa ctgcgaccag | 1800 |
| cgcacccgct agcagatcag gttaagatca taacacactc cggaagatca ggaaggtagc | 1860 |
| cggtcgaacc atacgacgct aaagtactga tgccagcagg aggtgccgta ccatggccag | 1920 |
| aattcctagc actgagtgag agcgccacgt tagtgtacaa cgaaagagag tttgtgaacc | 1980 |
| gcaaaactata ccacattgcc atgcatggcc ccgccaagaa tacagaagag gagcagtaca | 2040 |
| aggttacaaa ggcagagctt gcagaaacag agtacgtgtt tgacgtggac aagaagcggt | 2100 |
| gcgttaagaa ggaagaagcc tcaggtctgg tcctctcggg agaactgacc aaccctccct | 2160 |
| atcatgagct agctctggag ggactgaaga cccgacctgc ggtcccgtac aaggtcgaaa | 2220 |
| caataggagt gataggcaca ccggggctcg gcaagtcagc tattatcaag tcaactgtca | 2280 |
| cggcacgaga tcttgttacc agcggaaaga aagaaaattg tcgcgaaatt gaggccgacg | 2340 |
| tgctaagact gaggggtatg cagattacgt cgaagacagt agattcggtt atgctcaacg | 2400 |
| gatgccacaa agccgtagaa gtgctgtacg ttgacgaagc gttcgctgac cacgcaggag | 2460 |
| cactacttgc cttgattgct atcgtcaggg ccgcaagaa ggtagtacta tgcggagacc | 2520 |
| ccatgcaatg cggattcttc aacatgatgc aactaaaggc acatttcaat caccctgaaa | 2580 |
| aagacatatg caccaagaca ttctacaagt atatctcccg gcgttgacac cagccagtta | 2640 |
| cagctattgt atcgacactg cattacgatg gaaagatgaa aaccacgaac ccgtgcaaga | 2700 |
| agaacattga aatcgatatt acagggggcca caaagccgaa gccaggggat atcatcctga | 2760 |
| catgtttccg cgggtgggtt aagcaattgc aaatcgacta tcccggacat gaagtaatga | 2820 |
| cagccgctgc ctcacaaggc ctaaccagaa aaggagtgtg tgccgtccgg caaaaagtca | 2880 |
| atgaaaaccc actgtacgag atcacatcag agcatgtgaa cgtgttgctc acccgactg | 2940 |
| aggacaggct agtgtggaaa accttgacgg gcgacccatg gattaagcag cccactaaca | 3000 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|------------|------------|------------|-------------|-------------|------------|------|
| tacctaag | aaactttcag | gctactatag | aggactggga | agctgaacac | aagggaataa | 3060 |
| ttgctgcaat | aaacagcccc | actccccgtg | ccaatccgtt | cagctgcaag | accaacgttt | 3120 |
| gctgggcgaa | agcattggaa | ccgatactag | ccacggccgg | tatcgtactt | accggttgcc | 3180 |
| agtggagcga | actgttccca | cagtttgccg | atgacaaacc | acattcggcc | atttacgcct | 3240 |
| tagacgtaat | ttgcattaag | tttttcggca | tggacttgac | aagcggactg | ttttctaaac | 3300 |
| agagcatccc | actaacgtac | catcccgccg | attcagcgag | gccggtagct | cattgggaca | 3360 |
| acagcccagg | aaccgcgaag | tatgggtacg | atcacgccat | tgccgccgaa | ctctcccgtg | 3420 |
| gatttccggt | gttccagcta | gctgggaagg | gcacacaact | tgatttgagc | acggggagaa | 3480 |
| ccagagttat | ctctgcacag | cataacctgg | ttccgggtgaa | ccgcaatctt | cctcacgcct | 3540 |
| tagtccccga | gtacaaggag | aagcaaccgg | gcccgggtcaa | aaaattcttg | aaccagttca | 3600 |
| aacaccactc | agtacttggt | gtatcagagg | aaaaaattga | agctccccgt | aagagaatcg | 3660 |
| aatggatcgc | cccgattggc | atagccgggt | cagataagaa | ctacaacctg | gctttcgggt | 3720 |
| ttccgccgca | ggcacggtac | gacctgggtg | tcatcaacat | tggaactaaa | tacagaaacc | 3780 |
| accactttca | gcagtgcgaa | gaccatgcgg | cgaccttaaa | aaccctttcg | cgttcggccc | 3840 |
| tgaattgcct | taaccagga | ggcacccctc | tggtgaagtc | ctatggctac | gccgaccgca | 3900 |
| acagtgagga | cgtagtcacc | gctcttgcca | gaaagtttgt | caggggtgtc | gcagcgagac | 3960 |
| cagattgtgt | ctcaagcaat | acagaaatgt | acctgatttt | ccgacaacta | gacaacagcc | 4020 |
| gtacacggca | attcaccccg | caccatctga | attgcgtgat | ttcgtccgtg | tatgagggta | 4080 |
| caagagatgg | agttggagcc | gcgccgtcat | accgcaccaa | aaggggagaat | attgctgact | 4140 |
| gtcaagagga | agcagttgtc | aacgcagcca | atccgctggg | tagaccaggc | gaaggagtct | 4200 |
| gccgtgccat | ctataaacgt | tggccgacca | gttttaccga | ttcagccacg | gagacaggca | 4260 |
| ccgcaagaat | gactgtgtgc | ctaggaaaga | aagtgatcca | cgcggtcggc | cctgattttc | 4320 |
| ggaagcacc | agaagcagaa | gccttgaaat | tgctacaaaa | cgccctaccat | gcagtggcag | 4380 |
| acttagtaaa | tgaacataac | atcaagtctg | tcgccatttc | actgctatct | acaggcattt | 4440 |
| acgcagccgg | aaaagaccgc | cttgaagtat | cacttaactg | cttgacaacc | gcgctagaca | 4500 |
| gaactgacgc | ggacgtaacc | atctattgcc | tggataagaa | gtggaaggaa | agaatcgacg | 4560 |
| cggcactcca | acttaaggag | tctgtaacag | agctgaagga | tgaagatatg | gagatcgacg | 4620 |
| atgagttagt | atggatccat | ccagacagtt | gcttgaaggg | aagaaaggga | ttcagtacta | 4680 |
| caaaaggaaa | attgtattcg | tacttcgaag | gcaccaaatt | ccatcaagca | gcaaaagaca | 4740 |
| tggcgagat | aaaggtcctg | ttccctaatt | accaggaaa | taatgaacaa | ctgtgtgcct | 4800 |
| acatattggg | tgagaccatg | gaagcaatcc | gcgaaaagtg | cccggtcgac | cataaccctg | 4860 |
| cgtctagccc | gcccaaaacg | ttgccgtgcc | tttgcatgta | tgccatgacg | ccagaaaggg | 4920 |
| tccacagact | tagaagcaat | aacgtcaaag | aagttacagt | atgctcctcc | accccccttc | 4980 |
| ctaagcacia | aattaagaat | gttcagaagg | ttcagtcac | gaaagtagtc | ctgtttaatc | 5040 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|-------------|------------|------------|------------|------------|------|
| cgcacactcc | cgcattcggt | cccgcccgta | agtacataga | agtgccagaa | cagcctaccg | 5100 |
| ctcctcctgc | acaggccgag | gaggccccc | aagttgtagc | gacaccgtca | ccatctacag | 5160 |
| ctgataacac | ctcgcttgat | gtcacagaca | tctcactgga | tatggatgac | agtagcgaag | 5220 |
| gctcactttt | ttcgagcttt | agcggatcgg | acaactctat | tactagtatg | gacagttggt | 5280 |
| cgtcaggacc | tagttcacta | gagatagtag | accgaaggca | ggtggtggtg | gctgacgttc | 5340 |
| atgccgtcca | agagcctgcc | cctattccac | cgccaaggct | aaagaagatg | gcccgcctgg | 5400 |
| cagcggcaag | aaaagagccc | actccaccgg | caagcaatag | ctctgagtc | ctccacctct | 5460 |
| cttttggtgg | ggtatccatg | tccctcggat | caattttcga | cggagagacg | gcccgccagg | 5520 |
| cagcgggtaca | acccttgga | acaggcccca | cggatgtgcc | tatgtctttc | ggatcgtttt | 5580 |
| ccgacggaga | gattgatgag | ctgagccgca | gagtaactga | gtccgaacct | gtcctgtttg | 5640 |
| gatcatttga | accgggagaa | gtgaactcaa | ttatatcgtc | ccgatcagcc | gtatcttttc | 5700 |
| cactacgcaa | gcagagacgt | agacgcagga | gcaggaggac | tgaatactga | ctaaccgggg | 5760 |
| taggtgggta | catattttcg | acggacacag | gccctgggca | cttgcaaaag | aagtcggttc | 5820 |
| tgcagaacca | gcttacagaa | ccgaccttgg | agcgcgaatg | cctggaaaga | attcatgccc | 5880 |
| cggtgctcga | cacgtcgaaa | gaggaacaac | tcaaactcag | gtaccagatg | atgcccaccg | 5940 |
| aagccaacaa | aagtaggtac | cagtctcgta | aagtagaaaa | tcagaaagcc | ataaccactg | 6000 |
| agcgactact | gtcaggacta | cgactgtata | actctgccac | agatcagcca | gaatgctata | 6060 |
| agatcaccta | tccgaaacca | ttgtactcca | gtagcgtagc | ggcgaactac | tccgatccac | 6120 |
| agttcgctgt | agctgtctgt | aacaactatc | tgcatgagaa | ctatccgaca | gtagcatctt | 6180 |
| atcagattac | tgacgagtac | gatgcttact | tggatatggt | agacgggaca | gtcgcctgcc | 6240 |
| tggatactgc | aaccttctgc | cccgttaagc | ttagaagtta | cccgaacaaa | catgagtata | 6300 |
| gagccccgaa | tatccgcagt | gcggttccat | cagcgatgca | gaacacgcta | caaaatgtgc | 6360 |
| tcattgccgc | aactaaaaga | aattgcaacg | tcacgcagat | gcgtgaactg | ccaacactgg | 6420 |
| actcagcgac | attcaatgtc | gaatgctttc | gaaaatatgc | atgtaatgac | gagtattggg | 6480 |
| aggagttcgc | tcggaagcca | attaggatta | ccactgagtt | tgtcaccgca | tatgtagcta | 6540 |
| gactgaaagg | ccctaaggcc | gccgcactat | ttgcaaagac | gtataatttg | gtcccattgc | 6600 |
| aagaagtgcc | tatggataga | ttcgtcatgg | acatgaaaag | agacgtgaaa | gttacaccag | 6660 |
| gcacgaaaca | cacagaagaa | agaccgaaag | tacaagtgat | acaagccgca | gaacccttgg | 6720 |
| cgactgctta | cttatgcggg | attcaccggg | aattagtgcg | taggcttacg | gccgtcttgc | 6780 |
| ttccaaacat | tcacacgctt | tttgacatgt | cggcggagga | ttttgatgca | atcatagcag | 6840 |
| aacacttcaa | gcaaggcgac | ccggtacttg | agacggatat | cgcatcattc | gacaaaagcc | 6900 |
| aagacgacgc | tatggcggtta | accggtctga | tgatcttgga | ggacctgggt | gtggatcaac | 6960 |
| cactactcga | cttgatcgag | tgcgcctttg | gagaaatata | atccacccat | ctacctacgg | 7020 |
| gtactcgttt | taaattcggg | gcgatgatga | aatccggaat | gttcctcaca | ctttttgtca | 7080 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|------------|------------|-------------|-------------|------------|-------------|------|
| acacagtttt | gaatgtcggt | atcgccagca | gagtactaga | agagcggcgt | aaaacgtcca | 7140 |
| gatgtgcagc | gttcattggc | gacgacaaca | tcatacatgg | agtagtatct | gacaaagaaa | 7200 |
| tggctgagag | gtgcgccacc | tggctcaaca | tggagggttaa | gatcatcgac | gcagtcacg | 7260 |
| gtgagagacc | accttacttc | tgcggcggat | ttatcttgca | agattcggtt | actttccacag | 7320 |
| cgtgccgcgt | ggcggatccc | ctgaaaaggc | tgtttaagtt | gggtaaaccg | ctcccagccg | 7380 |
| acgacgagca | agacgaagac | agaagacgcg | ctctgctaga | tgaaacaaag | gcgtgggttta | 7440 |
| gagtaggtat | aacaggcact | ttagcagtgg | ccgtgacgac | ccggtatgag | gtagacaata | 7500 |
| ttacacctgt | cctactggca | ttgagaactt | ttgcccagag | caaaagagca | ttccaagcca | 7560 |
| tcagagggga | aataaagcat | ctctacgggt | gtcctaaata | gtcagcatag | tacatttcac | 7620 |
| ctgactaata | ctacaacacc | accaccatga | atagaggatt | ctttaacatg | ctcggccgcc | 7680 |
| gccccctccc | ggccccact | gccatgtgga | ggccgcggag | aaggaggcag | gcggccccga | 7740 |
| tgctgccccg | caacgggctg | gcttctcaaa | tccagcaact | gaccacagcc | gtcagtgtccc | 7800 |
| tagtcattgg | acaggcaact | agacctcaac | ccccacgtcc | acgcccgcc | ccgcgccaga | 7860 |
| agaagcaggc | gccaagcaa | ccaccgaagc | cgaagaaacc | aaaaacgcag | gagaagaaga | 7920 |
| agaagcaacc | tgcaaaaccc | aaaccggaa | agagacagcg | catggcactt | aagttggagg | 7980 |
| ccgacagatt | gttcgacgtc | aagaacgagg | acggagatgt | catcgggcac | gcactggcca | 8040 |
| tggaaggaaa | ggtaatgaaa | cctctgcacg | tgaaaggaa | catcgaccac | cctgtgctat | 8100 |
| caaagctcaa | atttaccaag | tcgtcagcat | acgacatgga | gttcgcacag | ttgccagtca | 8160 |
| acatgagaag | tgaggcattc | acctacacca | gtgaacaccc | cgaaggattc | tataactggc | 8220 |
| accacggagc | ggtgcagtat | agtggaggta | gattttaccat | ccctcgcgga | gtaggaggca | 8280 |
| gaggagacag | cggtcgtccg | atcatggata | actccggtcg | ggttgctcg | atagtcctcg | 8340 |
| gtggcgctga | tgaaggaaca | cgaactgccc | tttcggctgt | cacctggaat | agtaaagggg | 8400 |
| agacaattaa | gacgaccccc | gaagggacag | aagagtgggtc | cgagcacca | ctgggtcacgg | 8460 |
| caatgtgttt | gtcgggaaat | gtgagcttcc | catgcgaccg | cccggccaca | tgctataccc | 8520 |
| gcgaaccttc | cagagccctc | gacatccttg | aagagaacgt | gaaccatgag | gcctacgata | 8580 |
| ccctgtctaa | tgccatattg | cggtgcggat | cgctctggcag | aagcaaaaga | agcgtcattg | 8640 |
| acgactttac | cctgaccagc | ccctacttgg | gcacatgctc | gtactgccac | catactgtac | 8700 |
| cgtgcttcag | ccctgttaag | atcgagcagg | tctgggacga | agcggacgat | aacaccatac | 8760 |
| gcatacagac | ttccgcccag | tttgatac | accaaagcgg | agcagcaagc | gcaaacaagt | 8820 |
| accgctacat | gtcgcttaag | caggatcaca | ccgttaaaga | aggcaccatg | gatgacatca | 8880 |
| agattagcac | ctcaggaccg | tgtagaaggc | ttagctacaa | aggatacttt | ctcctcgcaa | 8940 |
| aatgccctcc | aggggacagc | gtaacgggtta | gcatagtgag | tagcaactca | gcaacgtcat | 9000 |
| gtacactggc | ccgcaagata | aaaccaaatt | tcgtgggacg | ggaaaaatat | gatctacctc | 9060 |
| ccgttcacgg | taaaaaatt | ccttgacacag | tgtacgaccg | tctgaaagaa | acaactgcag | 9120 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|-------------|------------|-------------|-------------|-------------|-------|
| gctacatcac | tatgcacagg | ccgagaccgc | acgcttatac | atcctacctg | gaagaatcat | 9180 |
| cagggaaagt | ttacgcaaag | ccgccatctg | ggaagaacat | tacgtatgag | tgcaagtgcg | 9240 |
| gcgactacaa | gaccggaacc | gtttcgaccc | gcaccgaaat | cactggttgc | accgccatca | 9300 |
| agcagtgcgt | cgcctataag | agcgacaaaa | cgaagtgggt | cttcaactca | ccggacttga | 9360 |
| tcagacatga | cgaccacacg | gccaagggga | aattgcattt | gcctttcaag | ttgatcccgga | 9420 |
| gtacctgcat | ggtcctgtt | gcccacgcgc | cgaatgtaat | acatggcttt | aaacacatca | 9480 |
| gcctccaatt | agatacagac | cacttgacat | tgtcaccac | caggagacta | ggggcaaacc | 9540 |
| cggaaccaac | cactgaatgg | atcgctggaa | agacggtcag | aaacttcacc | gtcgaccgag | 9600 |
| atggcctgga | atacatatgg | ggaaatcatg | agccagttag | ggcttatgcc | caagagttag | 9660 |
| caccaggaga | ccctcacgga | tggccacacg | aaatagtaca | gcattactac | catcgccatc | 9720 |
| ctgtgtacac | catcttagcc | gtcgcacacg | ctaccgtggc | gatgatgatt | ggcgtaactg | 9780 |
| ttgcagtgtt | atgtgcctgt | aaagcgcgcc | gtgagtgcct | gacgccatac | gccctggccc | 9840 |
| caaacgccgt | aatcccaact | tcgctggcac | tcttgtgctg | cgtagggtcg | gccaatgctg | 9900 |
| aaacgttcac | cgagaccatg | agttacttgt | ggtcgaacag | tcagccgttc | ttctgggtcc | 9960 |
| agttgtgcat | acctttggcc | gctttcatcg | ttctaattcg | ctgctgtctc | tgctgcctgc | 10020 |
| cttttttagt | ggttgccggc | gcctacctgg | cgaaggtaga | cgcctacgaa | catgcgacca | 10080 |
| ctgttccaaa | tgtgccacag | ataccgtata | aggcacttgt | tgaaggggca | gggtatgccc | 10140 |
| cgctcaattt | ggagatcact | gtcatgtcct | cggaggtttt | gccttccacc | aaccaagagt | 10200 |
| acattacctg | caaattcacc | actgtggtcc | cctcccaaaa | aatcaaatgc | tgcggtcctt | 10260 |
| tggaaatgtca | gccggccgct | catgcagact | atacctgcaa | ggtcttcgga | ggggtctacc | 10320 |
| cctttatgtg | gggaggagcg | caatgttttt | gcgacagtga | gaacagccag | atgagttagg | 10380 |
| cgtacgtcga | attgtcagca | gattgcgcgt | ctgaccacgc | gcaggcgatt | aagggtcaca | 10440 |
| ctgccgcgat | gaaagtagga | ctgcgtattg | tgtacgggaa | cactaccagt | ttccttagatg | 10500 |
| tgtacgtgaa | cggagtcaca | ccaggaacgt | ctaaagactt | gaaagtcata | gctggaccaa | 10560 |
| tttcagcatc | gtttacgcca | ttcgatcata | aggctgttat | ccatcgcggc | ctggtgtaca | 10620 |
| actatgactt | cccggaatat | ggagcgatga | aaccaggagc | gtttggagac | attcaagcta | 10680 |
| cctccttgac | tagcaaggat | ctcatcgcca | gcacagacat | taggctactc | aagccttccg | 10740 |
| ccaagaacgt | gcatgtcccg | tacacgcagg | cctcatcagg | atttgagatg | tggaaaaaca | 10800 |
| actcaggccg | cccactgcag | gaaaccgcac | ctttcggggtg | taagattgca | gtaaatccgc | 10860 |
| tccgagcggg | ggactgttca | tacgggaaca | ttcccatttc | tattgacatc | ccgaacgctg | 10920 |
| cctttatcag | gacatcagat | gcaccactgg | tctcaacagt | caaagtgtgaa | gtcagttagt | 10980 |
| gcacttattc | agcagacttc | ggcgggatgg | ccaccctgca | gtatgtatcc | gaccgcgaag | 11040 |
| gtcaatgccc | cgtacattcg | cattcgagca | cagcaactct | ccaagagtcg | acagtacatg | 11100 |
| tcctggagaa | aggagcgggtg | acagtacact | ttagcaccgc | gagtccacag | gcgaacttta | 11160 |

CI0042PCTseqlisting.ST25

tcgtatcgct gtgtgggaag aagacaacat gcaatgcaga atgtaaacca ccagctgacc 11220
 atatcgtagag caccgccgac aaaaatgacc aagaatttca agccgccatc tcaaaaacat 11280
 catggagttg gctgtttgcc cttttcgcg ggcctcgtc gctattaatt ataggactta 11340
 tgatttttgc ttgcagcatg atgctgacta gcacacgaag atgaccgcta cgccccaatg 11400
 atccgaccag caaaactcga tgtacttccg aggaactgat gtgcataatg catcaggctg 11460
 gtacattaga tccccgctta ccgcgggcaa tatagcaaca ctaaaaactc gatgtacttc 11520
 cgaggaagcg cagtgcataa tgctgcgag tgttgccaca taaccactat attaaccatt 11580
 tatctagcgg acgccaaaaa ctcaatgtat ttctgaggaa gcgtggtgca taatgccacg 11640
 cagcgtctgc ataactttta ttatttcttt tattaatcaa caaaattttg tttttaacat 11700
 ttc 11703

<210> 5
 <211> 10945
 <212> DNA
 <213> West Nile virus

<400> 5
 gaggattaac aacaattaac acagtgcgag ctgtttctta gcacgaagat ctcatgtct 60
 aagaaaccag gagggcccg caagagccgg gctgtcaata tgctaaaacg cggaatgccc 120
 cgctgtttgt ctttgattgg actgaagagg gctatgttga gcctgatcga cggcaagggg 180
 ccaatacgat ttgtgttggc tctcttggcg ttcttcaggt tcacagcaat tgctccgacc 240
 cgagcagtgc tggatcgatg gagaggtgtg aacaaacaaa cagcgatgaa acaccttctg 300
 agttttaaga aggaactagg gaccttgacc agtgcattca atcggcggag ctcaaaacaa 360
 aagaaaagag gaggaagac cggaattgca gtcattgatt gcctgatcgc cagcgtagga 420
 gcagttaccc tctctaactt ccaaggggaag gtgatgatga cggtaaatgc tactgacgtc 480
 acagatgtca tcacgattcc aacagctgct ggaaagaacc tatgcattgt cagagcaatg 540
 gatgtgggat acatgtgcga tgatactatc acctatgaat gccagtgct gtcggctggt 600
 aatgatccag aagacatcga ctgttggtgc acaaagtcag cagtctacgt caggtatgga 660
 agatgcacca agacacgcca ctcaagacgc agtcggaggt cactgacagt gcagacacac 720
 ggagaaagca ctctagcgaa caagaagggg gcttggtgag acagcaccaa ggccacaagg 780
 tatttggtaa aaacagaatc atggatcttg aggaaccctg gatatgccct ggtggcagcc 840
 gtcattggtt ggatgcttgg gagcaacacc atgcagagag ttgtgtttgt cgtgctattg 900
 cttttggtgg ccccgactta cagcttcaac tgccttgga tgagcaacag agacttcttg 960
 gaaggagtgt ctggagcaac atgggtggat ttggttctcg aaggcgatag ctgcgtgact 1020
 atcatgtcta aggacaagcc taccatcgat gtgaagatga tgaatatgga ggcggccaac 1080
 ctggcagagg tccgcagtta ttgctatttg gctaccgtca gcgatctctc caccaaagct 1140
 gcgtgcccga ccatggggga agcccacaat gacaaacgtg ctgacccagc ttttgtgtgc 1200
 agacaaggag tgggtggacag gggctggggc aacggctgcg gactatttgg caaaggaagc 1260

CI0042PCTseqlisting.ST25

| | |
|--|------|
| attgacacat gcgccaaatt tgcctgctct accaaggcaa taggaagaac catcttgaaa | 1320 |
| gagaatatca agtacgaagt ggccattttt gtccatggac caactactgt ggagtcgcac | 1380 |
| ggaaactact ccacacaggt tggagccact caggcagggg gattcagcat cactcctgcg | 1440 |
| gcgccttcat acacactaaa gcttggagaa tatggagagg tgacagtgga ctgtgaacca | 1500 |
| cggtcagggg ttgacaccaa tgcatactac gtgatgactg ttggaacaaa gacgttcttg | 1560 |
| gtccatcgtg agtggttcat ggacctcaac ctcccttggg gcagtgtctgg aagtactgtg | 1620 |
| tggaggaaca gagagacggt aatggagttt gaggaaccac acgccacgaa gcagtctgtg | 1680 |
| atagcattgg gctcacaaga gggagctctg catcaagctt tggctggagc cattcctgtg | 1740 |
| gaattttcaa gcaacactgt caagttgacg tcgggtcatt tgaagtgtag agtgaagatg | 1800 |
| gaaaaattgc agttgaaggg aacaacctat ggcgtctgtt caaaggcttt caagtttctt | 1860 |
| gggactcccg cagacacagg tcacggcact gtggtgttgg aattgcagta cactggcacg | 1920 |
| gatggacctt gcaaagttcc tatctcgtca gtggcttcat tgaacgacct aacgccagtg | 1980 |
| ggcagattgg tcaactgtcaa cccttttggt tcagtggcca cggccaacgc taaggctctg | 2040 |
| attgaattgg aaccaccctt tggagactca tacatagtgg tgggcagagg agaacaacag | 2100 |
| atcaatcacc attggcacia gtctggaagc agcattggca aagcctttac aaccaccctc | 2160 |
| aaaggagcgc agagactagc cgctctagga gacacagctt gggactttgg atcagttgga | 2220 |
| ggggtgttca cctcagttgg gaaggctgtc catcaagtgt tcggaggagc attccgctca | 2280 |
| ctgttcggag gcatgtcctg gataacgcaa ggattgtctg gggctctcct gttgtggatg | 2340 |
| ggcatcaatg ctcgtgatag gtccatagct ctcacgtttc tcgcagttgg aggagtctg | 2400 |
| ctcttcctct ccgtgaacgt gcacgctgac actgggtgtg ccataaacat caccgcaa | 2460 |
| gagctgagat gtggaagtgg agtgttcata cacaatgatg tggaggcttg gatggaccgg | 2520 |
| tacaagtatt accctgaaac gccacaaggc ctagccaaga tcattcagaa agctcataag | 2580 |
| gaaggagtgt gcggtctacg atcagtttcc agactggagc atcaaagtgt ggaagcagtg | 2640 |
| aaggacgagc tgaacactcc tttgaaggag aatggtgtgg acctagtgt cgtggttgag | 2700 |
| aaacaggagg gaatgtacaa gtcagcacct aaacgcctca ccgccaccac ggaaaaattg | 2760 |
| gaaattggct ggaaggcctg gggaaagagt attttatttg caccagaact cgccaacaac | 2820 |
| acctttgtgg ttgatggtcc ggagaccaag gaatgtccga ctcagaatcg cgcttggaat | 2880 |
| agcttagaag tggaggattt tggatttggg ctcaccagca ctcggatgtt cctgaaggtc | 2940 |
| agagaaggca acacaactga atgtgactcg aagatcattg gaacggctgt caagaacaac | 3000 |
| ttggcgatcc acagtgcctt gtcctatttg attgaaagca ggctcaatga tacgtggaag | 3060 |
| cttgaaaggg cagttctggg tgaagtcaaa tcattgtacgt ggcctgagac gcataccttg | 3120 |
| tggggcgatg gaatccttga gagtgacttg ataataccag tcacactggc gggaccacga | 3180 |
| agcaatcaca atcggagacc tgggtacaag acacaaaacc agggcccatg ggacgaaggc | 3240 |
| cgggtagaga ttgacttcga ttactgcccc ggaactacgg tcaccctgag tgagagctgc | 3300 |

CI0042PCTseqlisting.ST25

| | |
|--|------|
| ggacaccgtg gacctgccac tcgcaccacc acagagagcg gaaagttgat aacagattgg | 3360 |
| tgctgcagga gctgcacctt accaccactg cgctaccaa ctgacagcgg ctgttggtat | 3420 |
| ggatggaga tcagaccaca gagacatgat gaaaagaccc tcgtgcagtc acaagtgaat | 3480 |
| gcttataatg ctgatatgat tgaccctttt cagttgggcc ttctggtcgt gttcttggcc | 3540 |
| accaggagg tccttcgcaa gaggtggaca gccaatgca gcatgccagc tatactgatt | 3600 |
| gctctgctag tcctggtgtt tgggggcatt acttacactg atgtgttacg ctatgtcatc | 3660 |
| ttgtggggg cagctttcgc agaattcaat tcgggaggag acgtggtaca cttggcgctc | 3720 |
| atggcgacct tcaagataca accagtgttt atggtggcat cgtttctcaa agcgagatgg | 3780 |
| accaaccagg agaactttt gttgatgttg gcggctgttt tctttcaaat ggcttatcac | 3840 |
| gatgcccgcc aaattctgct ctgggagatc cctgatgtgt tgaattcact ggcggtagct | 3900 |
| tggatgatac tgagagccat aacattcaca acgacatcaa acgtggttgt tccgctgcta | 3960 |
| gccctgctaa caccgggct gagatgcttg aatctggatg tgtacaggat actgctgttg | 4020 |
| atggtcggaa taggcagctt gatcaggag aagaggagt cagctgcaa aaagaaagga | 4080 |
| gcaagtctgc tatgcttggc tctagcctca acaggacttt tcaaccccat gatccttgct | 4140 |
| gctggactga ttgcatgtga tccaaccgt aaacgcggat ggcccgaac tgaagtgatg | 4200 |
| acagctgtcg gcctaattgt tgccatcgtc ggagggttg cagagcttga cattgactcc | 4260 |
| atggccattc caatgactat cgcggggctc atgtttgctg ctttcgtgat ttctgggaaa | 4320 |
| tcaacagata tgtggattga gagaacggcg gacatttcct gggaaagtga tgcagaaatt | 4380 |
| acaggctcga gcgaaagagt tgatgtgcgg ctgatgatg atggaaactt ccagctcatg | 4440 |
| aatgatccag gagcaccttg gaagatatgg atgctcagaa tggctctgtc cgcgattagt | 4500 |
| gcgtacaccc cctgggcaat cttgccctca gtatgtggat ttgggataac tctccaatac | 4560 |
| acaaagagag gaggcgtgtt gtgggacact ccctcaccaa aggagtacaa aaagggggac | 4620 |
| accaccaccg gcgtttacag gatcatgact cgtgggctgc tcggcagtta tcaagcagga | 4680 |
| gcgggcgtga tggttgaagg tgttttccac accctttggc atacaacaaa aggagccgct | 4740 |
| ttgatgagcg gagagggccg cctggaccca tactggggca gtgtcaagga ggatcgactt | 4800 |
| tgttacggag gaccctggaa attgcagcac aagtggaacg ggcaggatga ggtgcagatg | 4860 |
| attgtggtgg aacctggcaa gaacgttaag aacgtccaga cgaaaccagg ggtgttcaaa | 4920 |
| acacctgaag gagaaatcgg ggccgtgact ttggacttcc cactggaac atcaggctca | 4980 |
| ccaatagtgg acaaaaacgg tgatgtgatt gggctttatg gcaatggagt cataatgcc | 5040 |
| aacggctcat acataagcgc gatagtgcag ggtgaaagga tggatgagcc aatcccagcc | 5100 |
| ggattcgaac ctgagatgct gaggaaaaaa cagatcactg tactggatct ccattcccggc | 5160 |
| gccggtaaaa caaggaggat tctgccacag atcatcaaa aggccataaa cagaagactg | 5220 |
| agaacagccg tgctagcacc aaccagggtt gtggctgctg agatggctga agcactgaga | 5280 |
| ggactgcca tccggtacca gacatccgca gtgccagag aacataatgg aatgagatt | 5340 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|------------|-------------|------------|------------|------------|-------------|------|
| gttgatgtca | tgtgtcatgc | taccctcacc | cacaggctga | tgtctcctca | caggggtgccg | 5400 |
| aactacaacc | tggtcgtgat | ggatgaggct | catttcaccg | accagctag | cattgcagca | 5460 |
| agaggttaca | tttccacaaa | ggtcgagcta | ggggaggcgg | cggcaatatt | catgacagcc | 5520 |
| acccaccag | gcacttcaga | tccattccca | gagtccaatt | caccaatttc | cgacttacag | 5580 |
| actgagatcc | cggatcgagc | ttggaactct | ggatacgaat | ggatcacaga | atacaccggg | 5640 |
| aagacggttt | ggtttgtgcc | tagtgtcaag | atggggaatg | agattgccct | ttgcctacaa | 5700 |
| cgtgctggaa | agaaagtagt | ccaattgaac | agaaagtcgt | acgagacgga | gtacccaaaa | 5760 |
| tgtaagaacg | atgattggga | ctttgttatc | acaacagaca | tatctgaaat | gggggctaac | 5820 |
| ttcaaggcga | gcaggggtgat | tgacagccgg | aagagtgtga | aaccaaccat | cataacagaa | 5880 |
| ggagaaggga | gagtgatcct | gggagaacca | tctgcagtga | cagcagctag | tgccgcccag | 5940 |
| agacgtggac | gtatcggtag | aaatccgtcg | caagttggtg | atgagtactg | ttatgggggg | 6000 |
| cacacgaatg | aagacgactc | gaacttcgcc | cattggactg | aggcacgaat | catgctggac | 6060 |
| aacatcaaca | tgccaaacgg | actgatcgct | caattctacc | aaccagagcg | tgagaaggta | 6120 |
| tataccatgg | atggggaata | ccggctcaga | ggagaagaga | gaaaaaactt | tctggaactg | 6180 |
| ttgaggactg | cagatctgcc | agtttggtcg | gcttacaagg | ttgcagcggc | tggagtgtca | 6240 |
| taccacgacc | ggaggtggtg | ctttgatggt | cctaggacaa | acacaatttt | agaagacaac | 6300 |
| aacgaagtgg | aagtcatcac | gaagcttggt | gaaaggaaga | ttctgaggcc | gcgctggatt | 6360 |
| gacgccaggg | tgtactcgga | tcaccaggca | ctaaaggcgt | tcaaggactt | cgcctcgggg | 6420 |
| aaacgttctc | agatagggct | cattgagggt | ctgggaaaga | tgcttgagca | cttcatgggg | 6480 |
| aagacatggg | aagcacttga | caccatgtac | gttgtggcca | ctgcagagaa | aggaggaaga | 6540 |
| gctcacagaa | tggccctgga | ggaactgcc | gatgctcttc | agacaattgc | cttgattgcc | 6600 |
| ttattgagtg | tgatgaccat | gggagtattc | ttctctctca | tgacgcggaa | gggcattgga | 6660 |
| aagataggtt | tgggaggcgc | tgtcttgga | gtagcgacct | tttctgttg | gatggctgaa | 6720 |
| gttccaggaa | cgaagatcgc | cggaatgttg | ctgctctccc | ttctcttgat | gattgtgcta | 6780 |
| attcctgagc | cagagaagca | acgttcgcag | acagacaacc | agctagccgt | gttcctgatt | 6840 |
| tgtgtcatga | cccttgtag | cgcagtggca | gccaacgaga | tgggttggt | agataagacc | 6900 |
| aagagtgaca | taagcagttt | gtttgggcaa | agaattgagg | tcaaggagaa | tttcagcatg | 6960 |
| ggagagtttc | ttctggactt | gaggccggca | acagcctggt | cactgtacgc | tgtgacaaca | 7020 |
| gcggtcctca | ctccactgct | aaagcatttg | atcacgtcag | attacatcaa | cacctcattg | 7080 |
| acctcaataa | acgttcaggc | aagtgcacta | ttcacactcg | cgcgaggctt | ccccttcgtc | 7140 |
| gatgttgagg | tgctggctct | cctgctagca | gccgatgct | ggggacaagt | cacctcacc | 7200 |
| gttacggtaa | cagcggcaac | actccttttt | tgccactatg | cctacatggt | tcccgtttgg | 7260 |
| caagctgagg | caatgcgctc | agcccagcgg | cggacagcgg | ccggaatcat | gaagaacgct | 7320 |
| gtagtgatg | gcacgtggc | cacggacgtc | ccagaattag | agcgcaccac | acccatcatg | 7380 |

CI0042PCTseqlisting.ST25

| | |
|---|------|
| cagaagaaag ttggacagat catgctgac | 7440 |
| ccgtctgtga agacagtacg agaagccgga attttgatca cggccgcagc ggtgacgctt | 7500 |
| tgggagaatg gagcaagctc tgtttgaac gcaacaactg ccatcggact ctgccacatc | 7560 |
| atgcgtgggg gttggtgtc atgtctatcc ataacatgga cactcataaa gaacatggaa | 7620 |
| aaaccagggc taaaagagg tggggcaaaa ggacgcacct tgggagagg tggaaagaa | 7680 |
| agactcaacc agatgacaaa agaagagttc actaggtacc gcaaagaggc catcatcgaa | 7740 |
| gtcgatcgct cagcggcaaa acacgccagg aaagaaggca atgtcactgg agggcatcca | 7800 |
| gtctctaggg gcacagcaaa actgagatgg ctggtcgaac ggaggtttct cgaaccggtc | 7860 |
| ggaaaagtga ttgaccttgg atgtggaaga ggcggttggg gttactatat ggcaacccaa | 7920 |
| aaaagagtcc aagaagtcag agggtagaca aaggcggtc ccggacatga agagcccaa | 7980 |
| ctagtgcaaa gttatggatg gaacattgtc accatgaaga gtggagtgga tgtgttctac | 8040 |
| agaccttctg agtggttga caccctcctt tgtgacatcg gagagtcctc gtcaagtgtc | 8100 |
| gaggttgaag agcataggac gattcgggtc cttgaaatgg ttgaggactg gctgcaccga | 8160 |
| gggccaaggg aattttgcgt gaagggtgc tgcccctaca tgccgaaagt catagagaag | 8220 |
| atggagctgc tccaacgccg gtatgggggg ggactggtca gaaaccact ctcacggaat | 8280 |
| tccacgcacg agatgtattg ggtgagtcga gcttcaggca atgtggtaca ttcagtgaat | 8340 |
| atgaccagcc aggtgctcct aggaagaatg gaaaaagga cctggaaggg accccaatac | 8400 |
| gaggaagatg taaacttggg aagtggaacc agggcggtgg gaaaaccct gctcaactca | 8460 |
| gacaccagta aaatcaagaa caggattgaa cgactcaggc gtgagtacag ttcgacgtgg | 8520 |
| caccacgatg agaaccaccc atatagaacc tggaactacc acggcagtta tgatgtgaag | 8580 |
| cccacaggct ccgccagttc gctggtcaat ggagtggta ggctcctctc aaaaccatgg | 8640 |
| gacaccatca cgaatgttac caccatggcc atgactgaca ctactcctt cgggcagcag | 8700 |
| cgagtgttca aagagaagg tggacacgaaa gctcctgaac cgccagaagg agtgaagtac | 8760 |
| gtgctcaacg agaccacca ctggttgtgg gcgtttttgg ccagagaaaa acgtcccaga | 8820 |
| atgtgctctc gagaggaatt cataagaaag gtcaacagca atgcagcttt gggtgccatg | 8880 |
| tttgaagagc agaatcaatg gaggagcgcc agagaagcag ttgaagatcc aaaattttgg | 8940 |
| gagatggtgg atgaggagcg cgaggcacat ctgcgggggg aatgtcacac ttgcatttac | 9000 |
| aacatgatgg gaaagagaga gaaaaaacc ggagagttcg gaaaggccaa gggagcaga | 9060 |
| gccatttggg tcatgtggct cgagctcgc tttctggagt tcgaggctct gggttttctc | 9120 |
| aatgaagacc actggcttgg aagaaagaac tcaggaggag gtgtcgagg cttgggcctc | 9180 |
| caaaaactgg gttacatcct gcgtgaagtt ggcacccggc ctgggggcaa gatctatgct | 9240 |
| gatgacacag ctggctggga caccgcac acgagagctg acttgaaaa tgaagctaag | 9300 |
| gtgcttgagc tgcttgatgg ggaacatcgg cgtcttgcca gggccatcat tgagctcacc | 9360 |
| tatcgtcaca aagttgtgaa agtgatgcgc ccggctgctg atggaagaac cgtcatggat | 9420 |

CI0042PCTseqlisting.ST25

```

gttatctcca gagaagatca gagggggagt ggacaagttg tcacctacgc cctaaacact 9480
ttcaccaacc tggccgtcca gctggtgagg atgatggaag gggaaaggagt gattggccca 9540
gatgatgtgg agaaactcac aaaagggaaa ggacccaaag tcaggacctg gctgtttgag 9600
aatggggaag aaagactcag ccgcatggct gtcagtggag atgactgtgt ggtaaagccc 9660
ctggacgata gctttgccac ctcgctccac ttcctcaatg ctatgtcaaa ggttcgcaaa 9720
gacatccaag agtggaaacc gtcaactgga tggatgatt ggcagcaggt tccattttgc 9780
tcaaaccatt tactgaatt gatcatgaa gatggaagaa cactggtggt tccatgccga 9840
ggacaggatg aattggtagg cagagctcgc atatctccag gggccggatg gaacgtccgc 9900
gacactgctt gtctggctaa gtcttatgcc cagatgtggc tgcttctgta cttccacaga 9960
agagacctgc ggctcatggc caacgccatt tgctccgctg tccctgtgaa ttgggtccct 10020
accggaagaa ccacgtggc catccatgca ggaggagagt ggatgacaac agaggacatg 10080
ttggaggtct ggaaccgtgt ttggatagag gagaatgaat ggatggaaga caaaacccca 10140
gtggagaaat ggagtacgt cccatattca ggaaaacgag aggacatctg gtgtggcagc 10200
ctgattggca caagagccc agccacgtgg gcagaaaaca tccagggtggc tatcaaccaa 10260
gtcagagcaa tcatcgga tgagaagtat gtggattaca tgagttcact aaagagatat 10320
gaagacacaa ctttggttga ggacacagta ctgtagatat ttaatcaatt gtaaatagac 10380
aatataagta tgcataaaag tgtagtttta tagtagtatt tagtggtgtt agtgtaaata 10440
gttaagaaaa ttttgaggag aaagtcaggc cggaagtgc ccgccaccg aagttgagta 10500
gacggtgctg cctgcgactc aaccccagga ggactgggtg acaaagccg cgaagtgatc 10560
catgtaagcc ctcaaacg tctcggaagg aggacccac atgttgtaac ttcaaagccc 10620
aatgtcagac cagctacgg cgtgctactc tgcggagagt gcagtctgcg atagtgtccc 10680
aggaggactg ggttaacaaa ggcaaaccaa cccccacgc ggccctagcc ccgtaatgg 10740
tgttaaccag ggcgaaagga ctagagggtta gaggagacc cgcggtttta agtgcacggc 10800
ccagcctggc tgaagctgta ggtcagggga aggactagag gttagtggag acccgtgcc 10860
acaaaacacc acaaaaaac agcatattga cacctgggat agactaggag atcttctgct 10920
ctgcacaacc agccacacgg cacag 10945

```

```

<210> 6
<211> 1542
<212> DNA
<213> Escherichia coli 16S Ribosomal RNA

```

```

<220>
<221> misc_feature
<222> (896)..(896)
<223> n is a, c, g, or t

```

```

<400> 6
aaattgaaga gtttgatcat ggctcagatt gaacgtggc ggcaggccta acacatgcaa 60

```

CI0042PCTseqlisting.ST25

gtcgaacggt aacaggaakc agcttgctga tttgctgacg agtggcggac gggtagagtaa 120
 tgtctgggaa actgcctgat ggagggggat aactactgga aacggtagct aataccgcat 180
 aacgtcgcaa gaccaaagag ggggaccttc gggcctcttg ccatcggaatg tgccagatg 240
 ggattagcta gtaggtgggg taaaggctca cctaggcgac gatccctagc tggctctgaga 300
 ggatgaccag ccacactgga actgagacac ggtccagact cctacgggag gcagcagtgg 360
 ggaatattgc acaatgggag caagcctgat gcagccatgc cgcgtgtatg aagaaggcct 420
 tcgggttgta aagtactttc agcggggagg aaggaggtaa agttaatacc tttgctcatt 480
 gacgttaccg gcagaagaag caccggctaa ctccgtgcca gcagccgagg taatacggag 540
 ggtgcaagcg ttaatcgga ttactgggag taaagcgac gcagggcgtt tgtaagtca 600
 gatgtgaaat ccccgggctc aacctgggaa ctgcatctga tactggcaag cttgagtctc 660
 gtagaggggg gtagaattcc aggtgtagcg gtgaaatgcg tagagatctg gaggaatacc 720
 ggtggcgaag gcggccccct ggacgaagac tgacgctcag gtgcgaaagc gtggggagca 780
 aacaggatta gataccctgg tagtccacgc cgtaaacgat gtcgacttgg aggttggtcc 840
 cttgaggcgt ggcttccgga gctaacgcgt taagtckacc gcctggggag tacgngcgca 900
 aggttaaaac tcaaataaat tgacgggggc cgcacaagc ggtggagcat gtggtttaat 960
 tcgatgcaac gcgaagaacc ttacctggtc ttgacatcca cagaactttc cagagatgga 1020
 ttggtgcctt cggaactgt gagacaggtg ctgcatggct gtcgtcagct cgtgttgtaga 1080
 aatgttgggt taagtccgc aacgagcgca acccttatcc tttgttgcca gcggtccggc 1140
 cggaactca aaggagactg ccagtataa actggaggaa ggtggggatg acgtcaagtc 1200
 atcatggccc ttacgaccag ggctacacac gtgctacaat ggcgcataca aagagaagcg 1260
 ayctcgag agcaagcgga cctcataaag tgcgtcgtag tccggattgg agtctgcaac 1320
 tcgactccat gaagtcggaa tcgctagtaa tcgtggatca gaatgccacg gtgaatacgt 1380
 tccccggcct tgtacacacc gcccgtcaca ccatgggagt ggggttgcaa agaagtaggt 1440
 agcttaacct tcgggagggc gcttaccact ttgtgattca tgactggggg gaagtcgtaa 1500
 caaggtaacc gtaggggaac ctgcggttgg atcacctcct ta 1542

<210> 7
 <211> 2905
 <212> DNA
 <213> Escherichia coli 23 S Ribosomal RNA

<400> 7
 ggttaagcga ctaagcgtag acggtggatg ccctggcagt cagaggcgat gaaggacgtg 60
 ctaatctgag ataagcgtag gtaaggatg atgaaccgtt ataaccggcg atttccgaat 120
 ggggaaaccc agtgtgtttc gacacactat cattaactga atccataggt taatgaggcg 180
 aaccggggga actgaaacat ctaagtaccc cgaggaaaag aaatcaaccg agattcccc 240
 agtagcggcg agcgaacggg gagcagccca gagcctgaat cagtgtgtgt gttagtggaa 300
 gcgtctggaa aggcgcgcga tacagggtga cagccccgta cacaaaaatg cacatgctgt 360

CI0042PCTseqlisting.ST25

| | |
|---|------|
| gagctcgatg agtagggcgg gacacgtggt atcctgtctg aatatggggg gaccatcctc | 420 |
| caaggctaaa tactcctgac tgaccgatag tgaaccagta ccgtgagggg aaggcgaaaa | 480 |
| gaaccccggc gaggggagtg aaaaagaacc tgaaccgtg tacgtacaag cagtgggagc | 540 |
| acgcttaggc gtgtgactgc gtaccttttg tataatgggt cagcgactta tattctgtag | 600 |
| caaggttaac cgaatagggg agccgaaggg aaaccgagtc ttaactgggc gttaagtgtc | 660 |
| agggtataga cccgaaaccc ggtgatctag ccatgggcag gttgaagggt gggtaacact | 720 |
| aactggagga ccgaaccgac taatgttgaa aaattagcgg atgacttgtg gctgggggtg | 780 |
| aaaggccaat caaacggga gatagctggt tctccccgaa agctatttag gtagcgctc | 840 |
| gtgaattcat ctccgggggt agagcactgt ttcggcaagg gggcatccc gacttacaa | 900 |
| cccgatgcaa actgcgaata ccggagaatg ttatcacggg agacacacgg cgggtgctaa | 960 |
| cgctccgtcg gaagaggga acaaccaga ccgccagcta aggtcccaa gtcattggtta | 1020 |
| agtgggaaac gatgtgggaa ggcccagaca gccaggatgt tggcttagaa gcagccatca | 1080 |
| tttaaagaaa gcgtaatagc tctgtggtcg agtcggcctg cgcggaagat gtaacggggc | 1140 |
| taaaccatgc accgaagctg cggcagcgac gcttatgctg tgttgggtag gggagcgtc | 1200 |
| tgtaagcctg cgaagggtgt ctgtgaggca tgctggagg atcagaagtg cgaatgctga | 1260 |
| cataagtaac gataaagcgg gtgaaaagcc cgctcgccgg aagaccaagg gttcctgtcc | 1320 |
| aacgttaatc ggggcagggt gagtcgaccc ctaaggcgag gccgaaaggc gtagtcgatg | 1380 |
| ggaaacaggt taatattcct gtacttggtg ttactgcgaa ggggggacgg agaaggctat | 1440 |
| gttggccggg cgacggttgt cccggtttaa gcgtgtaggc tggttttcca ggcaaaccg | 1500 |
| gaaaatcaag gctgaggcgt gatgacgagg cactacggtg ctgaagcaac aaatgccctg | 1560 |
| cttcaggaa aagcctctaa gcatcaggta acatcaaata gtaccccaa ccgacacagg | 1620 |
| tggtcaggta gagaatacca aggcgcttga gagaactcgg gtgaaggaa taggcaaaat | 1680 |
| ggtgccgtaa cttcgggaga aggcacgctg atatgtaggt gaagcgactt gctcgtggag | 1740 |
| ctgaaatcag tcgaagatac cagctggctg caactgttta ttaaaaacac agcactgtgc | 1800 |
| aaacacgaaa gtggacgtat acggtgtgac gcctgcccgg tgccggaagg ttaattgatg | 1860 |
| gggttagcgg taacgcgaag ctcttgatcg aagccccgt aaacggcggc cgtaactata | 1920 |
| acggctctaa ggtagcgaag ttccttgctg ggtaagttcc gacctgcacg aatggcgtaa | 1980 |
| tgatggccag gctgtctcca cccgagactc agtgaaattg aactcgctgt gaagatgcag | 2040 |
| tgtacccgcg gcaagacgga aagacccgt gaacctttac tatagcttga cactgaacat | 2100 |
| tgagccttga tgttaggat aggtgggagg ctttgaagtg tggacgccag tctgcatgga | 2160 |
| gccgaccttg aaataaccac ctttaattgt tgatgttcta acgttggccc gtaatccggg | 2220 |
| ttgcggacag tgtctggtgg gtagtttgac tggggcggtc tcctcctaaa gagtaacgga | 2280 |
| ggagcacgaa ggttggctaa tcctggctcg acatcaggag gttagtgcga tggcataagc | 2340 |
| cagcttgact gcgagcgtga cggcgcgagc aggtgcgaaa gcaggtcata gtgatccggt | 2400 |

CI0042PCTseqlisting.ST25

| | |
|--|------|
| ggttctgaat ggaagggcca tcgctcaacg gataaaaggt actccgggga taacaggctg | 2460 |
| ataccgcccc agagttcata tcgacggcgg tgtttggcac ctcgatgtcg gctcatcaca | 2520 |
| tcctggggct gaagtaggtc ccaaggggat ggctgttcgc catttaaagt ggtacgcgag | 2580 |
| ctgggttttag aacgtcgtga gacagttcgg tccctatctg ccgtgggcgc tggagaactg | 2640 |
| aggggggctg ctcctagtac gagaggaccg gagtggacgc atcactgggtg ttcgggttgt | 2700 |
| catgccaatg cactgcccgg tagctaaatg cggaagagat aagtgtgaa agcatctaag | 2760 |
| cacgaaactt gccccgagat gagttctccc tgacccttta agggctcctga aggaacgttg | 2820 |
| aagacgacga cgttgatagg ccgggtgtgt aagcgcagcg atgcgttgag ctaaccggta | 2880 |
| ctaatagaacc gtgaggctta acctt | 2905 |

<210> 8
 <211> 1798
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

| | |
|---|------|
| <400> 8 | |
| tatctggttg atcctgccag tagtcatatg cttgtctcaa agattaagcc atgcatgtct | 60 |
| aagtataagc aatttataca gtgaaactgc gaatggctca ttaaatacagt tatcgtttat | 120 |
| ttgatagttc ctttactaca tgggtataacc gtggtaattc tagagctaata acatgcttaa | 180 |
| aatctcgacc ctttgggaaga gatgtattta ttagataaaa aatcaatgtc ttcggactct | 240 |
| ttgatgattc ataataactt ttcgaatcgc atggccttgt gctggcgatg gttcattcaa | 300 |
| atcttctgcc tatcaacttt cgatggtagg atagtggcct accatggttt caacgggtaa | 360 |
| cggggaataa gggttcgatt ccggagaggg agcctgagaa acggctacca catccaagga | 420 |
| aggcagcagg cgcgcaaatt acccaatcct aattcaggga ggtagtgaca ataaataacg | 480 |
| atacagggcc cattcgggtc ttgtaattgg aatgagtaca atgtaaatac cttaacgagg | 540 |
| aacaattgga gggcaagtct ggtgccagca gccgcggtaa ttccagctcc aatagcgtat | 600 |
| attaaagtgt ttgcagttaa aaagctcgta gttgaacttt gggcccgggt ggccgggtccg | 660 |
| atcttttctgt gtactggatt tccaacgggg cctttccttc tggctaacct tgagtccttg | 720 |
| tggctcttgg cgaaccagga cttttacttt gaaaaaatta gagtgttcaa agcaggcgta | 780 |
| ttgctcgaat atattagcat ggaataatag aataggacgt ttggttctat tttgttggtt | 840 |
| tctaggacca tcgtaatgat taatagggac ggtcgggggc atcggtattc aattgtcgag | 900 |
| gtgaaattct tggatttatt gaagactaac tactgcgaaa gcatttgcca aggacgtttt | 960 |
| cattaatcaa gaacgaaagt taggggatcg aagatgatct ggtaccgtcg tagtcttaac | 1020 |
| cataaactat gccgactaga tcgggtggtg tttttttaat gaccactcg gtaccttacg | 1080 |
| agaaatcaaa gtctttgggt tctgggggga gtatggtcgc aaggctgaaa cttaaaggaa | 1140 |
| ttgacggaag ggcaccacta ggagtggagc ctgcggctaa tttgactcaa cacggggaaa | 1200 |
| ctcaccaggt ccagacacaa taaggattga cagattgaga gctctttctt gattttgtgg | 1260 |

CI0042PCTseqlisting.ST25

gtggtggtgc atggccgttt ctgagttggt ggagtgattt gtctgcttaa ttgcgataac 1320
 gaacgagacc ttaacctact aaatagtggt gctagcattt gctggttatt cacttcttag 1380
 aggactatc ggtttcaagc cgatggaagt ttgaggcaat aacaggctctg tgatgccctt 1440
 agaacgttct gggccgcacg cgcgctacac tgacggagcc agcaggtcta accttgccg 1500
 agaggctctg gtaatcttgt gaaactccgt cgtgctgggg atagagcatt gtaattattg 1560
 ctcttcaacg aggaattcct agtaagcgca agtcacgagc ttgcttgatg tacgtccctg 1620
 cctttgtac acaccgcccg tcgctagtac cgattgaatg gcttagtgag gcctcaggat 1680
 ctgcttagag aagggggcaa ctccatctca gagcggagaa tttggacaaa cttggtcatt 1740
 tagaggaact aaaagtcgta acaagggttc cgtaggtgaa cctgcggaag gatcatta 1798

<210> 9
 <211> 3911
 <212> DNA
 <213> Yeast 25S Ribosomal RNA

<400> 9
 aattccgtga tgggccttta ggttttacca actgcggcta atcttttttt atactgagcg 60
 tattggaacg ttatcgataa gaagagagcg tctaggcgaa caatgttctt aaagtttgac 120
 ctcaaactag gtaggagtag ccgctgaact taagcatatc aataagcggg ggaagagaaa 180
 ccaaccggat tgccttagta acggcgagtg aagcggcaaa agctcaaatt tgaaatctgg 240
 taccttcggt gcccgagttg taatttgag agggcaactt tggggccggt ccttgcttat 300
 gttccttgga acaggacgct atagaggggt agcatcccgt gtggcgagga gtgcggttct 360
 ttgtaaagtg ccttcgaaga gtcgagttgt ttgggaatgc agctctaagt gggtggtaaa 420
 ttccatctaa agctaaatat tggcgagaga ccgatagcga acaagtacag tgatggaaaag 480
 atgaaaagaa ctttgaaaag agagtgaaa agtacgtgaa attgttgaaa ggggaaggga 540
 tttgatcaga catggtgttt tgtgccctct gctccttggt ggtaggggaa tctcgcat 600
 cactgggcca gcatcagttt tgggtggcagg ataatccat aggaatgtag cttgcctcgg 660
 taagtattat agcctgtggg aatactgcca gctgggactg aggactgcga cgtaagtcaa 720
 ggatgctggc ataagtgtta tatgccgcc gtcttgaaac acggaccaag gagtctaacg 780
 tctatgcgag tgtttgggtg taaaacccat acgcgtaatg aaagtgaacg taggttgggg 840
 cctcgcaaga ggtgcacaa cgaccgatcc tgatgtcttc ggatggattt gagtaagagc 900
 atagctgttg ggacccgaaa gatggtgaac tatgcctgaa tagggtgaag ccagaggaaa 960
 ctctggtgga ggctcgtagc ggttctgacg tgcaaatcga tcgtcgaatt tgggtatagg 1020
 ggcgaaagac taatcgaacc atctagtagc tggttcctgc cgaagtttcc ctcaggatag 1080
 cagaagctcg tatcagtttt atgaggtaaa gcgaatgatt agaggttccg gggtcgaaat 1140
 gaccttgacc tattctcaaa ctttaaataat gtaagaagtc cttgttactt aattgaacgt 1200
 ggacatttga atgaagagct tttagtgggc ctttttggt aagcagaact ggcgatgcgg 1260
 gatgaaccga acgtagagtt aaggtgccgg aatacacgct catcagacac cacaaaaggt 1320

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|------------|------------|------------|------------|------------|------|
| gtagttcat | ctagacagcc | ggacggtggc | catggaagtc | ggaatccgct | aaggagtgtg | 1380 |
| taacaactca | ccggccgaat | gaactagccc | tgaaaatgga | tggcgctcaa | gcgtgttacc | 1440 |
| tatactctac | cgtcagggtt | gatatgatgc | cctgacgagt | aggcaggcgt | ggaggtcagt | 1500 |
| gacgaagcct | agaccgtaag | gtcgggtcga | acggcctcta | gtgcagatct | tgggtgtagt | 1560 |
| agcaaataatt | caaatagaaa | ctttgaagac | tgaagtgggg | aaaggttcca | cgtcaacagc | 1620 |
| agttggacgt | gggttagtcg | atcctaagag | atggggaagc | tccgtttcaa | aggcctgatt | 1680 |
| ttatgcaggc | caccatcgaa | agggaatccg | gtaagattcc | ggaacttggg | tatggattct | 1740 |
| tcacggtaac | gtaactgaat | gtggagacgt | cggcgcgagc | cctgggagga | gttatctttt | 1800 |
| cttcttaaca | gcttatcacc | ccggaattgg | tttatccgga | gatggggctt | tatggctgga | 1860 |
| agaggccagc | acctttgctg | gctccggtgc | gcttgtgacg | gcccgtgaaa | atccacagga | 1920 |
| aggaatagtt | ttcatgctag | gtcgtactga | taaccgcagc | aggtctcaa | ggtgaacagc | 1980 |
| ctctagttga | tagaataatg | tagataaggg | aagtcggcaa | aatagatccg | taacttcggg | 2040 |
| ataaggattg | gctctaaggg | tcgggtagtg | agggccttgg | tcagacgcag | cgggcgtgct | 2100 |
| tgtggactgc | ttggtggggc | ttgctctgct | aggcggacta | cttgcgtgcc | ttgtttaga | 2160 |
| cggccttggt | aggtctcttg | tagaccgtcg | cttgtacaaa | ttaacagatc | aacttagaac | 2220 |
| tggtacggac | aaggggaatc | tgactgtcta | attaaaacat | agcattgcga | tggtcagaaa | 2280 |
| gtgatgttga | cgcaatgtga | tttctgcccc | gtgctctgaa | tgtcaaagtg | aagaaattca | 2340 |
| accaagcgcg | agtaaacggc | gggagtaact | atgactctct | taaggtagcc | aatgcctcg | 2400 |
| tcacttaatt | agtgcgcgc | atgaatggat | taacgagatt | cccactgtcc | ctatctacta | 2460 |
| tctagcgaat | ccacagccaa | gggaacgggc | ttggcagaat | cagcggggaa | agaagaccct | 2520 |
| gttgagcttg | actctagttt | gacattgtga | agagacatag | agggtgtaga | ataagtggga | 2580 |
| gcttcggcgc | cagtgaataa | ccactacctt | tatagtttct | ttacttattc | aatgaagcgg | 2640 |
| agctggaatt | cattttccac | gttctagcat | tcaaggtccc | attcggggct | gatccggggt | 2700 |
| gaagacattg | tcaggtgggg | agtttggtcg | gggcggcaca | tctgttaaac | gataacgcag | 2760 |
| atgtcctaag | gggggctcat | ggagaacaga | aatctccagt | agaacaaaag | ggtaaagccc | 2820 |
| cttagtttga | tttcagtgtg | aatacaaac | attgaaagtg | tggcctatcg | atcctttagt | 2880 |
| ccctcggaat | ttgaggctag | agggtccaga | aaagttacca | cagggataac | tggcttgtgg | 2940 |
| cagtcaagcg | ttcatagcga | cattgctttt | tgattcttcg | atgtcggctc | ttcctatcat | 3000 |
| accgaagcag | aattcggtaa | gcgttggtat | gttcaccac | taatagggaa | catgagctgg | 3060 |
| gtttagaccg | tcgtgagaca | ggttagtttt | accctactga | tgaatgttac | cagcaatagt | 3120 |
| aattgaactt | agtacgagag | gaacagttca | ttcggataat | tggtttttgc | ggctgtctga | 3180 |
| tcaggcattg | ccgcgaagca | ccatccgctg | gattatggct | gaacgcctct | aagtcagaat | 3240 |
| ccatgctaga | acgcggtgat | ttctttgctc | cacacaatat | agatggatac | gaataaggcg | 3300 |
| tccttgtggc | gtcgtgtaac | catagcaggc | tagcaacggg | gcacttggcg | gaaaggcctt | 3360 |

CI0042PCTseq1isting.ST25

ggggtgcttg c tggcgaattg caatgtcatt ttgcgtgggg ataaatcatt tgtatacgac 3420
 ttagatgtac aacgggggtat tgtaagcggg agagtagcct tgttggttacg atctgctgag 3480
 attaacccctt tgttgcttga tttgtttttt atttctttct aagtgggtac tggcaggagc 3540
 cggggcctag ttttagagaga agtagactca acaagtctct ataaatttta tttgtcttaa 3600
 gaattctatg atccgggtaa aaacatgtat tgtatatac tattataata tacgatgagg 3660
 atgatagtgt gtaagagtgt accatttact aatgtatgta agttactatt tactatttgg 3720
 tctttttatt ttttattttt tttttttttt tcgttgcaaa gatgggttga aagagaaggg 3780
 ctttcacaaa gcttcccag cgtgaaagga tttgcccga cagtttgctt catggagcag 3840
 tttttccgc accatcagag cggcaaacat gagtgcttgt ataagtttag agaattgaga 3900
 aaagctcatt t 3911

<210> 10
 <211> 16569
 <212> DNA
 <213> Human mitochondrial DNA

<400> 10
 gatcacaggt ctatcacctt attaacctt cagggagct ctccatgcat ttggtatttt 60
 cgtctggggg gtatgcacgc gatagcattg cgagacgctg gagccggagc accctatgtc 120
 gcagtatctg tctttgattc ctgcctcatc ctattattta tcgcacctac gttcaatatt 180
 acaggcgaac atacttacta aagtgtgtta attaatatt gctttagtagga cataataata 240
 acaattgaat gtctgcacag ccactttcca cacagacatc ataacaaaaa atttccacca 300
 aacccccct ccccgctt tggccacagc acttaaacc atctctgcca aacccccaaa 360
 acaagaacc ctaaccag cctaaccaga tttcaaattt tatcttttgg cggtagcac 420
 ttttaacagt ccccccaa ctaacacatt attttcccct cccactcca tactactaat 480
 ctcatcaata caacccccgc ccactctacc cagcacacac acaccgctgc taacccata 540
 cccgaacca accaaacccc aaagacacc cccacagttt atgtagctta cctcctcaa 600
 gcaatacact gaaaatgtt agacgggctc acatcacccc ataaacaaat aggtttggtc 660
 ctagcctttc tattagctct tagtaagatt acacatgcaa gcatccccgt tccagttagt 720
 tcaccctcta aatcaccag atcaaaagg acaagcatca agcacgcagc aatgcagctc 780
 aaaacgctta gcctagccac acccccacgg gaaacagcag tgattaacct ttagcaataa 840
 acgaaagttt aactaagcta tactaacccc agggttggtc aatttcgtgc cagccaccgc 900
 ggtcacacga ttaaccaag tcaatagaag cggcgtaaa gagtgtttta gatcaccccc 960
 tccccataa agctaaaact cacctgagtt gtaaaaaact ccagttgaca caaaatagac 1020
 tacgaaagt gctttaacat atctgaacac acaatagcta agacccaaac tgggattaga 1080
 taccacctta tgcttagccc taaacctca cagttaaatc aacaaaactg ctcgccagaa 1140
 cactacgagc cacagcttaa aactcaaagg acctggcggg gcttcatatc cctctagagg 1200

CI0042PCTseqlisting.ST25

| | |
|---|------|
| agcctgttct gtaatcgata aaccccgatc aacctcacca cctcttgctc agcctatata | 1260 |
| ccgccatctt cagcaaacc tgaatgaaggc tacaagtaa gcgcaagtac ccacgtaaag | 1320 |
| acgttaggtc aaggtgtagc ccatgaggtg gcaagaaatg ggctacattt tctaccccag | 1380 |
| aaaactacga tagcccttat gaaacttaag ggtcgaagggt ggatttagca gtaaaactaag | 1440 |
| agtagagtgc ttagttgaac agggccctga agcgcgtaga caccgcccgt caccctctc | 1500 |
| aagtatactt caaaggacat ttaactaaaa cccctacgca tttatataga ggagacaagt | 1560 |
| cgtaacatgg taagtgtact ggaaagtgc cttggacgaa ccagagtgtg gcttaacaca | 1620 |
| aagcacccaa cttacactta ggagatttca acttaacttg accgctctga gctaaaccta | 1680 |
| gccccaaacc cactccacct tactaccaga caaccttagc caaacattt acccaaataa | 1740 |
| agtataggcg atagaaattg aaacctggcg caatagatat agtaccgcaa gggaaagatg | 1800 |
| aaaaattata accaagcata atatagcaag gactaacccc tataccttct gcataatgaa | 1860 |
| ttaactagaa ataactttgc aaggagagcc aaagctaaga ccccgaaac cagacgagct | 1920 |
| acctaagaac agctaaaaga gcacaccgt ctatgtagca aaatagtggg aagatttata | 1980 |
| ggtagaggcg acaaacctac cgagcctggt gatagctggt tgtccaagat agaacttag | 2040 |
| ttcaacttta aatttgccca cagaaccctc 'taaatcccct tgtaaattta actgttagtc | 2100 |
| caaagaggaa cagctctttg gacactagga aaaaaccttg tagagagagt aaaaaattta | 2160 |
| acacccatag taggcctaaa agcagccacc aattaagaaa gcgttcaagc tcaacacca | 2220 |
| ctacctaaaa aatcccaaac atataactga actcctcaca cccaattgga ccaatctatc | 2280 |
| accctataga agaactaatg ttagtataag taacatgaaa acattctcct ccgcataagc | 2340 |
| ctgcgtcaga ttaaaactt gaactgacaa ttaacagccc aatatctaca atcaaccaac | 2400 |
| aagtcattat taccctcact gtcaaccaa cacaggcatg ctcataagga aagggttaaaa | 2460 |
| aaagtaaaag gaactcggca aatcttacc cgctgttta ccaaaaacat cacctctagc | 2520 |
| atcaccagta ttagaggcac cgctgcccc gtgacacatg tttaacggcc gcggtaccct | 2580 |
| aaccgtgcaa aggtagcata atcacttggt ccttaaatag ggacctgtat gaatggctcc | 2640 |
| acgagggttc agctgtctct tacttttaac cagtgaattt gacctgcccg tgaagaggcg | 2700 |
| ggcataacac agcaagacga gaagacccta tggagcttta atttattaat gcaaacagta | 2760 |
| cctaacaac ccacaggtcc taaactacca aacctgcatt aaaaatttcg gttggggcga | 2820 |
| cctcggagca gaaccaacc tccgagcagt acatgctaag acttcaccag tcaaagcgaa | 2880 |
| ctactatact caattgatcc aataacttga ccaacggaac aagttaccct agggataaca | 2940 |
| gcgcaatcct attctagagt ccatatcaac aatagggttt acgacctcga tgttgatca | 3000 |
| ggacatcccg atggtgcagc cgctattaaa ggttcgtttg ttcaacgatt aaagtcctac | 3060 |
| gtgatctgag ttcagaccgg agtaatccag gtcggtttct atctaccttc aaattcctcc | 3120 |
| ctgtacgaaa ggacaagaga aataaggcct acttcacaaa gcgccttccc ccgtaaatga | 3180 |
| tatcatctca acttagtatt ataccacac ccaccaaga acagggtttg ttaagatggc | 3240 |

CI0042PCTseqlisting.ST25

| | | | | | | |
|-------------|------------|-------------|------------|-------------|------------|------|
| agagcccggg | aatcgcataa | aacttaaaac | tttacagtca | gaggttcaat | tcctcttctt | 3300 |
| aacaacatac | ccatggccaa | cctcctactc | ctcattgtac | ccatttcta | cgcaatggca | 3360 |
| ttcctaatac | ttaccgaacg | aaaaattcta | ggctatatac | aactacgcaa | aggccccaac | 3420 |
| gttgtaggcc | cctacgggct | actacaaccc | ttcgtgacg | ccataaaact | cttcaccaa | 3480 |
| gagcccctaa | aacccgccac | atctaccatc | accctctaca | tcaccgcccc | gaccttagct | 3540 |
| ctcaccatcg | ctcttctact | atgaaccccc | ctccccatac | ccaaccccc | ggtcaacctc | 3600 |
| aacctaggcc | tcctatttat | tctagccacc | tctagcctag | ccgtttactc | aatcctctga | 3660 |
| tcagggtag | catcaaactc | aaactacgcc | ctgatcggcg | caactgcgagc | agtagcccaa | 3720 |
| acaatctcat | atgaagtcac | cctagccatc | attctactat | caacattact | aataagtggc | 3780 |
| tcctttaacc | ttccaccctt | tatcacaaca | caagaacacc | tctgattact | cctgccatca | 3840 |
| tgacccttgg | ccataatatg | atttatctcc | acactagcag | agaccaaccg | aacccccctc | 3900 |
| gaccttgccg | aaggggagtc | cgaactagtc | tcaggcttca | acatcgaata | cgccgcaggc | 3960 |
| cccttcgccc | tattcttcat | agccgaatac | acaacatta | ttataataaa | caccctcacc | 4020 |
| actacaatct | tcctaggaac | aacatatgac | gcactctccc | ctgaactcta | cacaacatat | 4080 |
| tttgtcacca | agaccctact | tctaacctcc | ctgttcttat | gaattcgaac | agcatacccc | 4140 |
| cgattccgct | acgaccaact | catacacctc | ctatgaaaaa | acttcctacc | actcaccta | 4200 |
| gcattactta | tatgatattg | ctccataacc | attacaatct | ccagcattcc | ccctcaaacc | 4260 |
| taagaaatat | gtctgataaa | agagttactt | tgatagagta | aataatagga | gcttaaacc | 4320 |
| cccttatttct | aggactatga | gaatcgaacc | catccctgag | aatccaaaat | tctccgtgcc | 4380 |
| acctatcaca | ccccatccta | aagtaaggct | agctaaataa | gctatcgggc | ccataccccg | 4440 |
| aaaatgttgg | ttataccctt | cccgtactaa | ttaatccctt | ggcccaacc | gtcatctact | 4500 |
| ctaccatctt | tgcaggcaca | ctcatcacag | cgctaagctc | gcactgattt | tttacctgag | 4560 |
| taggcctaga | aataaacatg | ctagctttta | ttccagttct | aaccaaaaaa | ataaacctc | 4620 |
| gttcacaga | agctgccatc | aagtatttcc | tcacgcaagc | aaccgcatcc | ataatcctt | 4680 |
| taatagctat | cctcttcaac | aatatactct | ccggacaatg | aaccataacc | aatactacca | 4740 |
| atcaatactc | atcattaata | atcataatag | ctatagcaat | aaaactagga | atagccccct | 4800 |
| ttcacttctg | agtcccagag | gttaccacaag | gcacccctct | gacatccggc | ctgcttcttc | 4860 |
| tcacatgaca | aaaactagcc | cccatctcaa | tcataataca | aatctctccc | tcactaaacg | 4920 |
| taagccttct | cctcactctc | tcaatcttat | ccatcatagc | aggcagttga | ggtggattaa | 4980 |
| accagacca | gctacgcaa | atcttagcat | actcctcaat | taccacata | ggatgaataa | 5040 |
| tagcagttct | accgtacaac | cctaacataa | ccattcttaa | tttaactatt | tatattatcc | 5100 |
| taactactac | cgcattccta | ctactcaact | taaactccag | caccacgacc | ctactactat | 5160 |
| ctcgcacctg | aaacaagcta | acatgactaa | cacccttaat | tccatccacc | ctcctctccc | 5220 |
| taggaggcct | gcccccgcta | accggctttt | tgcccaaatg | ggccattatc | gaagaattca | 5280 |

CI0042PCTseqlisting.ST25

| | |
|--|------|
| caaaaaacaa tagcctcatc atccccacca tcatagccac catcacccctc cttaacctct | 5340 |
| acttctacct acgcctaatac tactccacct caatcacact actccccata tctaacaacg | 5400 |
| taaaaataaa atgacagttt gaacatacaa aaccacccc attcctcccc acactcatcg | 5460 |
| cccttaccac gctactccta cctatctccc cttttatact aataatctta tagaaattta | 5520 |
| ggttaaatac agaccaagag ctttcaaagc cctcagtaag ttgcaatact taattttctgt | 5580 |
| aacagctaag gactgcaaaa cccactctg catcaactga acgcaaatca gccactttaa | 5640 |
| ttaagctaag cccttactag accaatggga cttaaaccga caaacactta gttaacagct | 5700 |
| aagcaccta atcaactggc ttcaatctac ttctccgcc gccgggaaaa aaggcgggag | 5760 |
| aagccccggc aggtttgaag ctgcttcttc gaatttgcaa ttcaatatga aaatcacctc | 5820 |
| ggagctggtg aaaagaggcc taaccctgt ctttagattt acagtccaat gcttcaactca | 5880 |
| gccattttac ctacccccca ctgatgttcg ccgaccgttg actattctct acaaaccaca | 5940 |
| aagacattgg aacactatac ctattattcg gcgcatgagc tggagtccta ggcacagctc | 6000 |
| taagcctcct tattcgagcc gagctgggccc agccaggcaa cttcttaggt aacgaccaca | 6060 |
| tctacaacgt tatcgtcaca gcccattgc tttgtaataat cttcttcata gtaataccca | 6120 |
| tcataatcgg aggcctttggc aactgactag ttcccctaata aatcgggtgcc cccgatatgg | 6180 |
| cgtttccccc cataaacaac ataagcttct gactcttacc tccctctctc ctactcctgc | 6240 |
| tcgcatctgc tatagtggag gccggagcag gaacagggtg aacagtctac cctcccttag | 6300 |
| cagggaacta ctcccaccct ggagcctccg tagacctaac catcttctcc ttacacctag | 6360 |
| cagggtgtctc ctctatctta ggggccatca atttcatcac aacaattatc aatataaaac | 6420 |
| cccctgccat aaccacaatac caaacgcccc tcttcgtctg atccgtccta atcacagcag | 6480 |
| tcctacttct cctatctctc ccagtcctag ctgctggcat cactatacta ctaacagacc | 6540 |
| gcaacctcaa caccaccttc ttcgaccccc ccggaggagg agacccatt ctataccaac | 6600 |
| acctattctg atttttcggt caccctgaag tttatatctt tatcctacca ggcttcggaa | 6660 |
| taatctccca tattgtaact tactactccg gaaaaaaga accatttggg tacataggta | 6720 |
| tggcttgagc tatgatatca attggcttcc tagggtttat cgtgtgagca caccatatat | 6780 |
| ttacagtagg aatagacgta gacacacgag catatttcac ctccgctacc ataatcatcg | 6840 |
| ctatccccac cggcgtcaaa gtatttagct gactcgccac actccacgga agcaatatga | 6900 |
| aatgatctgc tgcagtgtc tgagccctag gattcatctt tcttttcacc gtagggtggcc | 6960 |
| tgactggcat tgtattagca aactcatcac tagacatcgt actacacgac acgtactacg | 7020 |
| ttgtagccca cttccactat gtcctatcaa taggagctgt atttgccatc ataggaggct | 7080 |
| tcattcactg atttccccta ttctcaggct acaccctaga ccaaacctac gccaaaatcc | 7140 |
| atttcactat catattcatc ggcgtaaatc taactttctt cccacaacac tttctcggcc | 7200 |
| tatccggaat gccccgacgt tactcggact accccgatgc atacaccaca tgaaacatcc | 7260 |
| tatcatctgt aggcctcattc atttctctaa cagcagtaat attaataatt ttcattgattt | 7320 |

CI0042PCTseqlisting.ST25

gagaagcctt cgcttcgaag cgaaaagtcc taatagtaga agaaccctcc ataaacctgg 7380
 agtgactata tggatgcccc ccaccctacc acacattcga agaaccgta tacataaaat 7440
 ctagacaaaa aaggaaggaa tcgaaccccc caaagctggt ttcaagccaa ccccatggcc 7500
 tccatgactt tttcaaaaag gtattagaaa aaccatttca taactttgtc aaagttaaata 7560
 tataggctaa atcctatata tcttaatggc acatgcagcg caagtaggtc tacaagacgc 7620
 tacttcccc atcatagaag agcttatcac ctttcatgat cagccctca taatcatttt 7680
 ccttatctgc ttcctagtcc tgtatgccct tttcctaaca ctcaacaata aactaactaa 7740
 tactaacatc tcagacgctc aggaaataga aaccgtctga actatcctgc ccgccatcat 7800
 cctagtctc atcgccctcc catccctacg catcctttac ataacagacg aggtcaacga 7860
 tccctccctt accatcaaat caattggcca ccaatggtac tgaacctacg agtacaccga 7920
 ctacggcgga ctaatcttca actcctacat acttccccca ttattcctag aaccaggcga 7980
 cctgcgactc cttagcgttg acaatcgagt agtactcccg attgaagccc ccattcgtat 8040
 aataattaca tcacaagacg tcttgactc atgagctgtc cccacattag gcttaaaaac 8100
 agatgcaatt cccggacgct taaaccaaac cactttcacc gctacacgac cgggggtata 8160
 ctacgggtcaa tgctctgaaa tctgtggagc aaaccacagt ttcattgccc tcgtcctaga 8220
 attaattccc ctaaaaatct ttgaaatagg gcccgatttt accctatagc accccctcta 8280
 cccctctag agccactgt aaagctaact tagcattaac cttttaagtt aaagattaag 8340
 agaaccaaca cctctttaca gtgaaatgcc ccaactaaat actaccgtat ggcccaccat 8400
 aattaccccc atactcctta cactattcct catcacccaa ctaaaaatat taaacacaaa 8460
 ctaccaccta cctccctcac caaagcccat aaaataaaaa aattataaca aaccctgaga 8520
 accaaaatga acgaaaatct gttcgcttca ttcatgccc ccacaatcct aggcctaccc 8580
 gccgcagtac tgatcattct atttccccct ctattgatcc ccacctcaa atatctcatc 8640
 aacaaccgac taatcaccac ccaacaatga ctaatcaaac taacctcaa acaaatgata 8700
 accatacaca aactaaagg acgaacctga tctcttatac tagtatcctt aatcattttt 8760
 attgccacaa ctaacctcct cggactcctg cctcactcat ttacaccaac cacccaacta 8820
 tctataaacc tagccatggc catcccccta tgagcgggca cagtgattat aggctttcgc 8880
 tctaagatta aaaatgccct agccacttc ttaccacaag gcacacctac accccttatc 8940
 cccatactag ttattatcga aaccatcagc ctactcatc aaccaatagc cctggccgta 9000
 cgctaaccg ctaacattac tgcaggccac ctactcatgc acctaattgg aagcgccacc 9060
 ctagcaatat caaccattaa ccttccctct acacttatca tcttcacaat tctaattcta 9120
 ctgactatcc tagaaatcgc tgtcgcttca atccaagcct acgttttcac acttctagta 9180
 agcctctacc tgcacgacaa cacataatga cccaccaatc acatgcctat catatagtaa 9240
 aaccagccc atgacccta acaggggccc tctcagccct cctaatagacc tccggcctag 9300
 ccatgtgatt tcacttccac tccataacgc tcctcatact aggcctacta accaacacac 9360

CI0042PCTseqlisting.ST25

| | |
|--|-------|
| taaccatata ccaatgatgg cgcgatgtaa cagcagaaaag cacataccaa ggccaccaca | 9420 |
| caccacctgt ccaaaaaggc cttcgatacg ggataatcct atttattacc tcagaagttt | 9480 |
| ttttcttcgc aggatttttc tgagcctttt accactccag cctagcccct acccccacat | 9540 |
| taggagggca ctggcccccacacagggcatca ccccgctaaa tcccctagaa gtcccactcc | 9600 |
| taaacacatc cgtattactc gcatcaggag tatcaatcac ctgagctcac catagtctaa | 9660 |
| tagaaaacaa ccgaaaccaa ataattcaag cactgcttat tacaatttta ctgggtctct | 9720 |
| attttaccct cctacaagcc tcagagtact tcgagtctcc cttcaccatt tccgacggca | 9780 |
| tctacggctc aacatttttt gtagccacag gcttcacgg acttcacgtc attattggct | 9840 |
| caactttcct cactatctgc ttcacccgcc aactaatatt tcactttaca tccaaacatc | 9900 |
| actttggctt cgaagccgcc gcctgatact ggcattttgt agatgtgggt tgactatttc | 9960 |
| tgtatgtctc catctattga tgagggtctt actcttttag tataaatagt accgttaact | 10020 |
| tccaattaac tagttttgac aacattcaaa aaagagtaat aaacttcgcc ttaattttaa | 10080 |
| taatcaacac cctcctagcc ttactactaa taattattac attttgacta ccacaactca | 10140 |
| acggctacat agaaaaatcc accccttacg agtgcggtt cgaccctata tccccgccc | 10200 |
| gcgtcccttt ctccataaaa ttcttcttag tagctattac cttcttatta tttgatctag | 10260 |
| aaattgccct ctttttacct ctaccatgag ccctacaaac aactaacctg ccactaatag | 10320 |
| ttatgtcatc cctcttatta atcatcatcc tagccctaag tctggcctat gagtgactac | 10380 |
| aaaaaggatt agactgaacc gaattggtat atagttttaa caaacgaat gatttcgact | 10440 |
| cattaaatta tgataatcat atttacaaa tgcccctcat ttacataaat attatactag | 10500 |
| catttaccat ctactttcta ggaatactag tatatcgctc acacctcata tcctccctac | 10560 |
| tatgcctaga aggaataata ctatcgctgt tcattatagc tactctcata accctcaaca | 10620 |
| cccactccct cttagccaat attgtgccta ttgccatact agtctttgcc gcctgcgaag | 10680 |
| cagcgggtggg cctagcccta ctagtctcaa tctccaacac atatggccta gactacgtac | 10740 |
| ataacctaaa cctactccaa tgctaaaact aatcgctcca acaattatat tactaccact | 10800 |
| gacatgactt tccaaaaaac acataatttg aatcaacaca accacccaca gcctaattat | 10860 |
| tagcatcatc cctctactat tttttaacca aatcaacaac aacctattta gctgttcccc | 10920 |
| aaccttttcc tccgaccccc taacaacccc cctcctaata ctaactacct gactcctacc | 10980 |
| cctcacaatc atggcaagcc aacgccactt atccagtga cactatcac gaaaaaaact | 11040 |
| ctacctctct atactaatct ccctacaaat ctcttaatt ataacattca cagccacaga | 11100 |
| actaatcata ttttatatct tcttcgaaac cacacttatc cccaccttgg ctatcatcac | 11160 |
| ccgatgaggc aaccagccag aacgcctgaa cgagggcaca tacttcctat tctacaccct | 11220 |
| agtaggctcc cttcccctac tcatcgact aatttacact cacaacaccc taggctcact | 11280 |
| aaacattcta ctactcactc tctactgccc agaactatca aactcctgag ccaataactt | 11340 |
| aatatgacta gcttacacaa tagcttttat agtaaagata cctctttacg gactccactt | 11400 |

CI0042PCTseqlisting.ST25

atgactccct aaagcccatg tcgaagcccc catcgctggg tcaatagtag ttgccgcagt 11460
 actcttaaaa ctaggcggct atggtataat acgcctcaca ctatttctca accccctgac 11520
 aaaacacata gcctaccctt tccttgtagt atccctatga ggcataatta taacaagctc 11580
 catctgccta cgacaaacag acctaaaatc gctcattgca tactcttcaa tcagccacat 11640
 agccctcgta gtaacagcca ttctcatcca aacccttgga agcttcaccg gcgcagtcac 11700
 tctcataatc gccacgggct ttacatcctc attactattc tgccatagcaa actcaaacta 11760
 cgaacgcact cacagtcgca tcataatcct ctctcaagga cttcaaactc tactcccact 11820
 aatagctttt tgatgacttc tagcaagcct cgctaacctc gccttaccct ccactattaa 11880
 cctactggga gaactctctg tgctagtaac caggttctcc tgatcaaata tcactctcct 11940
 acttacagga ctcaacatac tagtcacagc cctatactcc ctctacatat ttaccacaac 12000
 acaatggggc tcactcacc accacattaa caacataaaa ccctcattca cagagaaaaa 12060
 caccctcatg ttcatacacc tatcccccat tctcctccta tccctcaacc ccgacatcat 12120
 taccgggttt tcctcttgta aatatagttt aacaaaaaca tcagattgtg aatctgacaa 12180
 cagaggctta cgaccctta tttaccgaga aagctcaca gaactgctaa ctcatgcccc 12240
 catgtctaac aacatggctt tctcaacttt taaaggataa cagctatcca ttggtcttag 12300
 gccccaaaaa ttttggtgca actccaaata aaagtaataa ccatgcacac tactataacc 12360
 accctaacc tgacttcctt aattcccccc atccttacca ccctcgtaa ccctaacaaa 12420
 aaaaactcat accccatta tgtaaaatcc attgtcgcat ccacctttat tatcagtctc 12480
 ttccccaca caatattcat gtgcctagac caagaagtta ttatctcgaa ctgacactga 12540
 gccacaacc aaacaacca gctctcccta agcttcaaac tagactactt ctccataata 12600
 ttcatccctg tagcattggt cggtacatgg tccatcatag aattctcact gtgatata 12660
 aactcagacc caaacattaa tcagttcttc aaatatctac tcatcttctt aattaccata 12720
 ctaatcttag ttaccgctaa caacctattc caactgttca tcggctgaga gggcgtagga 12780
 attatctctt tcttgctcat cagttgatga tacgcccag cagatgcaa cacagcagcc 12840
 attcaagcaa tcctatacaa ccgtatcggc gatatcgggt tcatcctcgc cttagcatga 12900
 tttatcctac actccaactc atgagaccca caacaaatag cccttctaaa cgctaattcca 12960
 agcctcacc cactactagg ctcctccta gcagcagcag gcaaatcagc ccaattaggt 13020
 ctccaccct gactccctc agcatagaa ggccccacc cagtctcagc cctactccac 13080
 tcaagcacta tagttgtagc aggaatcttc ttactcatcc gcttccacc cctagcagaa 13140
 aatagccac taatccaaac tctaacta tgcttaggcg ctatcaccac tctgttcgca 13200
 gcagtctgcg cccttacaca aaatgacatc aaaaaaatcg tagccttctc cacttcaagt 13260
 caactaggac tcataatagt tacaatcggc atcaaccaac cacacctagc attcctgcac 13320
 atctgtacc acgccttctt caaagccata ctatttatgt gctccgggtc catcatccac 13380
 aaccttaaca atgaacaaga tattcgaaaa ataggaggac tactcaaac catacctctc 13440

CI0042PCTseqlisting.ST25

acttcaacct ccctcaccat tggcagccta gcattagcag gaataccttt cctcacaggt 13500
 ttctactcca aagaccacat catcgaaacc gcaaacatat catacacaaa cgctgagcc 13560
 ctatctatta ctctcatcgc tacctccctg acaagcgct atagcactcg aataattctt 13620
 ctcaccctaa caggtaacc tcgcttcccc acccttacta acattaacga aaataacccc 13680
 accctactaa accccattaa acgcctggca gccggaagcc tattcgagg atttctcatt 13740
 actaacaaca tttccccgc atcccccttc caaacaacaa tccccctta cctaaaactc 13800
 acagccctcg ctgtcacttt ctaggactt ctaacagccc tagacctcaa ctacctaacc 13860
 aacaaactta aaataaaatc cccactatgc acattttatt tctccaacat actcggattc 13920
 taccctagca tcacacaccg cacaatcccc tatctaggcc ttcttacgag ccaaaacctg 13980
 cccctactcc tcctagacct aacctgacta gaaaagctat tacctaaaac aatttcacag 14040
 caccaaatct ccacctccat catcacctca acccaaaaag gcataattaa actttacttc 14100
 ctctctttct tcttccact catcctaacc ctactcctaa tcacataacc tattcccccg 14160
 agcaatctca attacaatat atacaccaac aaacaatggt caaccagtaa ctactactaa 14220
 tcaacgccc taatcatata aagccccgc accaatagga tcctcccgaa tcaaccctga 14280
 cccctctct tcataaatta ttcagcttcc tacactatta aagtttacca caaccaccac 14340
 cccatcatac tctttcacc acagaccaa tcctacctcc atcgctaacc ccactaaaac 14400
 actcaccaag acctcaacc ctgaccccca tgctcagga tactcctcaa tagccatcgc 14460
 tgtagtatat ccaaagacaa ccattctcc ccctaataa attaaaaaaa ctattaaacc 14520
 catataacct ccccaaaaat tcagaataat aacacacccg accacaccgc taacaatcaa 14580
 tactaaaccc ccataaatag gagaaggctt agaagaaaac cccacaaacc ccattactaa 14640
 acccacactc aacagaaaca aagcatacat cattattctc gcacggacta caaccacgac 14700
 caatgatatg aaaaaccatc gttgtatttc aactacaaga acaccaatga cccaatacg 14760
 caaaattaa cccctaataa aattaattaa cactcattc atcgacctcc ccacccatc 14820
 caacatctcc gcatgatgaa acttcggctc actccttggc gcctgcctga tcctccaaat 14880
 caccacagga ctattcctag ccattgacta ctcaccagac gcctcaaccg cttttctatc 14940
 aatcgccac atcactcgag acgtaaatta tggctgaatc atccgctacc ttcacgcaa 15000
 tggcgctca atattcttta tctgcctctt cctacacatc gggcgaggcc tatattacgg 15060
 atcatttctc tactcagaaa cctgaaacat cggcattatc ctctgcttg caactatagc 15120
 aacagccttc ataggctatg tcctcccggtg aggccaaata tcattctgag gggccacagt 15180
 aattacaac ttactatccg ccattccata cattgggaca gacctagt cctaatctg 15240
 aggaggctac tcagtagaca gtccaccct cacacgattc ttacctttc acttcatctt 15300
 gcccttcatt attgcagccc tagcaacact ccacctccta ttcttgacg aaacgggatc 15360
 aaacaacccc cttagaatca cctccattc cgataaaatc accttcacc cttactacac 15420
 aatcaaagac gccctcggct tacttctctt ccttctctcc ttaatgacat taacactatt 15480

CI0042PCTseqlisting.ST25

```

ctcaccagac ctcctaggcg acccagacaa ttatacccta gccaacccct taaacacccc 15540
tccccacatc aagcccgaat gatatttcct attcgccctac acaattctcc gatccgtccc 15600
taacaaacta ggaggcgctc ttgccctatt actatccatc ctcacccctag caataatccc 15660
catcctccat atatccaaac aacaaagcat aatatttcgc ccactaagcc aatcacttta 15720
ttgactccta gccgcagacc tcctcattct aacctgaatc ggaggacaac cagtaagcta 15780
cccttttacc atcattggac aagtagcatc cgtactatac ttcacaacaa tcctaatacct 15840
aataccaact atctccctaa ttgaaaacaa aataactcaa tgggcctgtc cttgtagtat 15900
aaactaatac accagtcttg taaaccggag atgaaaacct tttccaagg acaaatacaga 15960
gaaaaagtct ttaactccac cattagcacc caaagctaag attctaattt aaactattct 16020
ctgttctttc atggggaagc agatttgggt accaccaag tattgactca cccatcaaca 16080
accgctatgt atttcgtaca ttactgccag ccaccatgaa tattgtacgg taccataaat 16140
acttgaccac ctgtagtaca taaaaaccca atccacatca aaacccctc cccatgctta 16200
caagcaagta cagcaatcaa cctcaacta tcacacatca actgcaactc caaagccacc 16260
cctcaccac taggatacca acaaacctac ccacccttaa cagtacatag tacataaagc 16320
catttaccgt acatagcaca ttacagtcaa atcccttctc gtcccatgg atgaccccc 16380
tcagataggg gtcccttgac caccatcctc cgtgaaatca atatcccgca caagagtgt 16440
actctcctcg ctccggggcc ataacacttg ggggtagcta aagtgaactg tatccgacat 16500
ctggttccta cttcaggggc ataaagccta aatagccac acgttcccct taaataagac 16560
atcacgatg                                     16569

```

<210> 11
<211> 17
<212> DNA
<213> B19 virus

<400> 11
tggtctggga tgaaggt 17

<210> 12
<211> 23
<212> DNA
<213> B19 virus

<400> 12
ccattttagg cgggcaaccc acc 23

<210> 13
<211> 21
<212> DNA
<213> B19 virus

<400> 13
tggaagtgt gctgtgcctg g 21

<210> 14
<211> 23

CI0042PCTseqlisting.ST25

<212> DNA
<213> B19 virus

<400> 14
agaatcattt gtcggaagct cag 23

<210> 15
<211> 20
<212> DNA
<213> B19 virus

<400> 15
acaagcctgg gcaagttagc 20

<210> 16
<211> 21
<212> DNA
<213> B19 virus

<400> 16
acaatgccag tggaaaggag g 21

<210> 17
<211> 23
<212> DNA
<213> B19 virus

<400> 17
ctttaacaca tgaagaccat gca 23

<210> 18
<211> 28
<212> DNA
<213> B19 virus

<400> 18
cctctcaaaa cactagaata tccttacg 28

<210> 19
<211> 27
<212> DNA
<213> B19 virus

<400> 19
cagccatacc accactggga cacagat 27

<210> 20
<211> 23
<212> DNA
<213> B19 virus

<400> 20
aatgccattt ctcatgggtca gac 23

<210> 21
<211> 20
<212> DNA
<213> B19 virus

<400> 21
gcaaaacaac accacaggca 20

CI0042PCTseqlisting.ST25

<210> 22
 <211> 22
 <212> DNA
 <213> Hepatitis B virus

 <400> 22
 caacctccaa tcactcacca ac 22

 <210> 23
 <211> 24
 <212> DNA
 <213> Hepatitis B virus

 <400> 23
 cctccaattt gtcctgggta tcgc 24

 <210> 24
 <211> 23
 <212> DNA
 <213> Hepatitis B virus

 <400> 24
 gtgtctgcgg cgttttatca tat 23

 <210> 25
 <211> 23
 <212> DNA
 <213> Hepatitis B virus

 <400> 25
 tgtttggtt tcagctatat gga 23

 <210> 26
 <211> 22
 <212> DNA
 <213> Hepatitis B virus

 <400> 26
 aattgtgggt cttttgggct tt 22

 <210> 27
 <211> 18
 <212> DNA
 <213> Hepatitis B virus

 <400> 27
 ttctccgtct gccgttcc 18

 <210> 28
 <211> 21
 <212> DNA
 <213> Hepatitis B virus

 <400> 28
 accagcacca tgcaactttt t 21

 <210> 29
 <211> 25
 <212> DNA

CI0042PCTseqlisting.ST25

<213> Hepatitis B virus
 <400> 29
 tttttcacct ctgcctaatac atctc 25

 <210> 30
 <211> 16
 <212> DNA
 <213> Hepatitis B virus
 <400> 30
 tcccactggt caagcc 16

 <210> 31
 <211> 21
 <212> DNA
 <213> Hepatitis B virus
 <400> 31
 acctcaccat accgcactca g 21

 <210> 32
 <211> 25
 <212> DNA
 <213> Hepatitis B virus
 <400> 32
 gggaattgat gactctagct acctg 25

 <210> 33
 <211> 22
 <212> DNA
 <213> Hepatitis B virus
 <400> 33
 atcctgaatg gcaaactcct tc 22

 <210> 34
 <211> 20
 <212> DNA
 <213> Hepatitis B virus
 <400> 34
 gagaaaccac acgtagcgca 20

 <210> 35
 <211> 16
 <212> DNA
 <213> Hepatitis B virus
 <400> 35
 caccctcca cacggc 16

 <210> 36
 <211> 27
 <212> DNA
 <213> Hepatitis B virus
 <400> 36
 atctctccac ctctaagaga cagtcac 27

CI0042PCTseqlisting.ST25

<210> 37
 <211> 22
 <212> DNA
 <213> Porcine parvovirus

 <400> 37
 cctcacaaaa cggcaagtac tg 22

 <210> 38
 <211> 36
 <212> DNA
 <213> Porcine parvovirus

 <400> 38
 acctaagtcc aagtgactgc tactgggttca tacagc 36

 <210> 39
 <211> 29
 <212> DNA
 <213> Porcine parvovirus

 <400> 39
 gttaataatg caatgcaaag tacctctaa 29

 <210> 40
 <211> 25
 <212> DNA
 <213> Porcine parvovirus

 <400> 40
 aaagacaact gaaagagagc atgga 25

 <210> 41
 <211> 30
 <212> DNA
 <213> Porcine parvovirus

 <400> 41
 ttcagctca gattctggct tcatgacaaa 30

 <210> 42
 <211> 22
 <212> DNA
 <213> Porcine parvovirus

 <400> 42
 caattctatt tcatgggcca gc 22

 <210> 43
 <211> 17
 <212> DNA
 <213> Porcine parvovirus

 <400> 43
 cgtggagcga gccaaaca 17

 <210> 44
 <211> 28
 <212> DNA
 <213> Porcine parvovirus

CI0042PCTseqlisting.ST25

| | | |
|---------------------------------|--|----|
| <400> 44 | | |
| ctgcacttaa ctccaacacc gccagatt | | 28 |
| <210> 45 | | |
| <211> 24 | | |
| <212> DNA | | |
| <213> Porcine parvovirus | | |
| <400> 45 | | |
| gcaatacgga caccaagtcc aact | | 24 |
| <210> 46 | | |
| <211> 23 | | |
| <212> DNA | | |
| <213> Porcine parvovirus | | |
| <400> 46 | | |
| gaggttaagaa gatcgccgag aaa | | 23 |
| <210> 47 | | |
| <211> 28 | | |
| <212> DNA | | |
| <213> Porcine parvovirus | | |
| <400> 47 | | |
| aacctcacca ccaaccaaaa tatataat | | 28 |
| <210> 48 | | |
| <211> 29 | | |
| <212> DNA | | |
| <213> Porcine parvovirus | | |
| <400> 48 | | |
| actactaact gaacctacca cagaaggag | | 29 |
| <210> 49 | | |
| <211> 26 | | |
| <212> DNA | | |
| <213> Porcine parvovirus | | |
| <400> 49 | | |
| cttttacctt cagatccaat aggagg | | 26 |
| <210> 50 | | |
| <211> 18 | | |
| <212> DNA | | |
| <213> Sindbis virus | | |
| <400> 50 | | |
| gcgtgcggac cctgtact | | 18 |
| <210> 51 | | |
| <211> 25 | | |
| <212> DNA | | |
| <213> Sindbis virus | | |
| <400> 51 | | |
| attggcttcg acaccacca gttca | | 25 |

CI0042PCTseqlisting.ST25

<210> 52
 <211> 20
 <212> DNA
 <213> Sindbis virus

 <400> 52
 ttctcggcta tggcaggttc 20

 <210> 53
 <211> 25
 <212> DNA
 <213> Sindbis virus

 <400> 53
 gtttatttct ccgtaggatc gacac 25

 <210> 54
 <211> 20
 <212> DNA
 <213> Sindbis virus

 <400> 54
 aaaactgctg caggtctcgg 20

 <210> 55
 <211> 22
 <212> DNA
 <213> Sindbis virus

 <400> 55
 gaaatcgata ttacaggggc ca 22

 <210> 56
 <211> 21
 <212> DNA
 <213> Sindbis virus

 <400> 56
 gcattaagtt tttcggcatg g 21

 <210> 57
 <211> 18
 <212> DNA
 <213> Sindbis virus

 <400> 57
 cgattggcat agccggtg 18

 <210> 58
 <211> 21
 <212> DNA
 <213> West Nile virus

 <400> 58
 tcagcgatct ctccaccaaa g 21

 <210> 59
 <211> 22
 <212> DNA
 <213> West Nile virus

CI0042PCTseqlisting.ST25

<400> 59
 tgcccgaacca tgggggaagc cc 22

<210> 60
 <211> 21
 <212> DNA
 <213> West Nile virus

<400> 60
 caatgacaaa cgtgctgacc c 21

<210> 61
 <211> 20
 <212> DNA
 <213> West Nile virus

<400> 61
 gctagtcctg gtgtttgggg 20

<210> 62
 <211> 24
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA

<400> 62
 agagtttgat catggctcag attg 24

<210> 63
 <211> 24
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA

<400> 63
 ctggcggcag gcctaacaca tgca 24

<210> 64
 <211> 21
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA

<400> 64
 aataccgcat aacgtcgcaa g 21

<210> 65
 <211> 20
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA

<400> 65
 gatgcaacgc gaagaacctt 20

<210> 66
 <211> 23
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA

<400> 66
 gactgggggtg aagtcgtaac aag 23

<210> 67

CI0042PCTseqlisting.ST25

<211> 24
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA
 <400> 67
 gtaacaaggt aaccgtaggg gaac 24

<210> 68
 <211> 19
 <212> DNA
 <213> Escherichia coli 16S Ribosomal RNA
 <400> 68
 gcggttgat cacctcctt 19

<210> 69
 <211> 22
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 69
 ccgatagtga accagtaccg tg 22

<210> 70
 <211> 25
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 70
 atgttgaaaa attagcggat gactt 25

<210> 71
 <211> 18
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 71
 gcactgtttc ggcaaggg 18

<210> 72
 <211> 18
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 72
 gccggaagac caagggtt 18

<210> 73
 <211> 24
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 73
 ggccgtaact ataacggtcc taag 24

<210> 74
 <211> 26
 <212> DNA
 <213> Escherichia coli 23S Ribosomal RNA
 <400> 74

CI0042PCTseqlisting.ST25

gataagtgct gaaagcatct aagcac 26

<210> 75
 <211> 25
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 75
 ctgccagtag tcatatgctt gtctc 25

<210> 76
 <211> 36
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 76
 tacagtgaag ctgcgaatgg ctcattaaat cagtta 36

<210> 77
 <211> 27
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 77
 taatacatgc ttaaaatctc gaccctt 27

<210> 78
 <211> 24
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 78
 gtcttcggac tctttgatga ttca 24

<210> 79
 <211> 17
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 79
 gcagccgcgg taattcc 17

<210> 80
 <211> 24
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 80
 gctgaaactt aaaggaattg acgg 24

<210> 81
 <211> 23
 <212> DNA
 <213> Yeast (*S. cerevisiae*)

<400> 81
 tggaagtttg aggcaataac agg 23

<210> 82
 <211> 19

CI0042PCTseqlisting.ST25

<212> DNA
 <213> Yeast (*S. cerevisiae*)
 <400> 82
 tgaacctgcg gaaggatca 19

<210> 83
 <211> 24
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 83
 aagcatatca ataagcggag gaaa 24

<210> 84
 <211> 20
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 84
 ctctggtgga ggctcgtagc 20

<210> 85
 <211> 18
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 85
 aatggatggc gctcaagc 18

<210> 86
 <211> 25
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 86
 tgaaaatcca caggaaggaa tagtt 25

<210> 87
 <211> 21
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 87
 ctaagggtcg ggtagtgagg g 21

<210> 88
 <211> 20
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 88
 agaaattcaa ccaagcgga 20

<210> 89
 <211> 19
 <212> DNA
 <213> Yeast 25S Ribosomal RNA
 <400> 89
 atgtcatttt gcgtgggga 19

CI0042PCTseqlisting.ST25

<210> 90
<211> 23
<212> DNA
<213> Human mitochondrial DNA

<400> 90
gttcaccctc taaatcacca cga 23

<210> 91
<211> 23
<212> DNA
<213> Human mitochondrial DNA

<400> 91
caagcacgca gcaatgcagc tca 23

<210> 92
<211> 26
<212> DNA
<213> Human mitochondrial DNA

<400> 92
ggaaacagca gtgattaacc tttagc 26

<210> 93
<211> 28
<212> DNA
<213> Human mitochondrial DNA

<400> 93
gactacgaaa gtggctttaa catatctg 28

<210> 94
<211> 24
<212> DNA
<213> Human mitochondrial DNA

<400> 94
tagagtgctt agttgaacag ggcc 24

<210> 95
<211> 24
<212> DNA
<213> Human mitochondrial DNA

<400> 95
taggcgatag aaattgaaac ctgg 24

<210> 96
<211> 21
<212> DNA
<213> Human mitochondrial DNA

<400> 96
tttgtaaga tggcagagcc c 21

<210> 97
<211> 20
<212> DNA

CI004227Tseqlisting.ST25

<213> Human mitochondrial DNA
<400> 97
agaatcgaac ccatccctga 20

<210> 98
<211> 19
<212> DNA
<213> Human mitochondrial DNA
<400> 98
tttcaccgta ggtggcctg 19

<210> 99
<211> 20
<212> DNA
<213> Human mitochondrial DNA
<400> 99
aatcgctgtc gccttaatcc 20

<210> 100
<211> 19
<212> DNA
<213> Escherichia coli 23S Ribosomal RNA
<400> 100
tcctacggga ggcagcagt 19

<210> 101
<211> 23
<212> DNA
<213> Escherichia coli 23S Ribosomal RNA
<400> 101
cgtattaccg cggctgctgg cac 23